Biotechnology Biotechnologie

BIOTECH

LCMS-IT-TOF™

The LCMS-IT-TOF mass spectrometer combines a quadrupole ion trap and a time of flight mass analyzer (TOF), resulting in a system which provides MSⁿ information with high mass accuracy and mass resolution for LC detection.

It delivers a wealth of qualitative information from a sample in a single experiment, thus eliminating multiple analyses and avoiding splitting of samples between multiple instruments. The outstanding performance of the LCMS-IT-TOF is based on the temperature controlled TOF detector which eliminates the need for regular mass calibration over a long period of time. Short duty cycles and consistent data quality of the LCMS-IT-TOF allow extensive fragmentation of the various eluting components of a single HPLC separation. In addition to the structural information gathered from related MS/MS and MSⁿ experiments, the empirical formulae for each compound or fragment can be predicted.

Thanks to the auto tuning function and the easy-to-use LCMSsolution software, this instrument enables scientists to generate outstanding data within a very short time.

With multiple ionization sources, it will be a critical tool for the analysis of small as well as larger molecules in target compound identification and characterization, particularly in the fields of impurity analysis, metabolomics, proteomics, metabolite identification and biomarker discovery.

- high mass accuracy (less than 5 ppm) and mass resolution (greater than 10,000) for all MS modes (MS, MS², MS³; MSⁿ)
- rapid polarity switching to support metabolism and biomarker discovery programs
- high precursor ion selection with precursor ion resolution of the ion trap of 1,000 at m/z 1,000
- high sensitivity of detection
- large variety of ionization sources (ESI, APCI, APPI and nano spray)
- additional software packages used in metabolomics, proteomics as well as in formula prediction are available to support data evaluation, e.g. Formula Predictor software)

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

Biotechnologie

LCMS-2010EV

The LCMS-2010EV single stage quadrupole mass analyzer has evolved into the most flexible MS detection system for LC currently available.

It delivers unprecedented performance for a complete range of analytical challenges from integrated column switching at capillary flow rates to mass directed auto-purification at preparative flow rates.

- maximized sensitivity through optimized ion optics
- high speed mass range scanning electronics (6,000 amu/s)
- supports ESI, APCI and APPI source
- unique APCI support for low flow rate applications (25 $\mu\text{L}/\text{min}$)
- supported by PsiPort™ web based software for both analytical and preparative data review
- in mass directed auto-purification the fraction triggering can be optimized using both slope and threshold parameters in addition to AND/OR logic utilizing other HPLC detector signals

- supports analog output for triggering an external fraction collector in polarity switching mode

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AXIMA™ MALDI-TOF MS Systems

The AXIMA range of MALDI-TOF MS systems delivers uncompromised data quality for biochemical, chemical and life science researchers.

The open architecture integrates with other enabling technologies from Shimadzu Biotech such as chemical inkjet printing (CHIP-1000) for tissue imaging and LC-MALDI (AccuSpot), resulting in the most flexible MALDI technology platform currently available. The AXIMA MALDI-TOF MS systems create a new benchmark in mass accuracy and resolution with high sensitivity.

AXIMA-Assurance

Linear MALDI-TOF mass spectrometer for mass detection without the need for structural information. The instrument is ideally suited for high throughput QA/QC applications in the areas of oligonucleotides, protein detection, identification of microorganisms, polymer analysis and microsatellite DNA characterization.

AXIMA-Confidence

MALDI-TOF mass spectrometer with curved field reflectron. The unique patented curved field reflectron provides in seconds the most efficient seamless PSD peptide fragmentation data for protein identification.

Automatic seamless PSD can be combined with automated peptide mass fingerprint experiments for further identification or confirmation.

AXIMA-Performance

MALDI-TOF mass spectrometer for high energy CID MS/MS with high performance, patented monoPULSE™ ion gate providing market-leading MS/MS

precursor ion selection resolution. Combined with high resolution MS spectra, more information can be extracted even from mixtures and neighboring isotopic envelopes.

The unique combination of the novaCFR™ reflectron and the 20-K™ high energy collision cell means that you see all the MS/MS fragment ions produced in the instrument.

The zoomIMMO™ feature allows rapid zooming in on diagnostic immonium ions at the low mass end of the spectrum.

High resolution and accuracy are consistent across the MS/MS mass range with no stitching of MS/MS fragment and precursor ion spectra.

No grids are used in the optics, thereby ensuring high ion transmission.


Sensitivity is improved for reliable MS/MS results, guaranteeing confident identifications without complications.

AXIMA-QIT

Hybrid MALDI-TOF mass spectrometer combining quadrupole ion trap technology with the simplicity of a time of flight MALDI system.

The design enables researchers to use the power of MSⁿ and the accuracy and resolution of a TOF detection. The AXIMA-QIT provides researchers with a unique tool set for the characterization of protein sequences, proteoglycans and carbohydrates.

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Biotechnologie

CHIP-1000 Chemical Inkjet Printer Chemischer Inkjet Printer

The Chemical Inkjet Printer is a unique piezoelectric liquid delivery system designed for exploring new challenges in micro-scale chemistry. Using nanotechnology fluidics the CHIP can deliver approx. 100 picoliter droplets of bio-chemicals and solvents to a diverse range of applications in biomedical and diagnostic analysis including on-tissue experiments (profiling/imaging), micro-scale on membrane protein identification and nanowestern analysis.

- delivers picoliter volumes of reagents to specific target locations (typically 100 – 400 μm in diameter) using 4 piezoelectric printing devices
- surface targets can be mounted on a conventional 96 well plate format (including PVDF membranes for PMF analysis and nanowesterns)
- focus on specific surface targets opens new horizons in working with critical biological samples (e.g. brain tissue for tissue profiling/ imaging)

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Biotechnology


Biotechnologie

- supports the delivery of multiple chemistries to a specific target location, increasing the level of protein information (sequence coverage) without repeating sample separation
- integrated with Shimadzu Biotech AXIMA MALDI mass spectrometers for rapid protein identification from tissue section as well as PVDF membranes

The CHIP-1000 technology was developed by Shimadzu and Proteome Systems Limited.

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AccuSpot LC-MALDI

This high capacity LC spotting device was designed for online collection of fractions eluting from Capillary/ Nano LC. It combines the identification performance and automation of a MALDI-MS system with the separation performance and flexibility of an LC system. Integrated with 1D or 2D LC systems, the AccuSpot increases significantly the potential to automate challenging proteomic studies.

- supports micro-fractionation on AXIMA MALDI-MS targets, for 96 or 384 well plates; 9 micro titer plates can be loaded at a time (using either 96 well or 384 well format)
- multiple plates can be spotted continuously, nine plates without stopping
- the flow rate of the co-delivered MALDI matrix solution can be adjusted from 0.1 $\mu\text{L}/\text{min}$ to 50 $\mu\text{L}/\text{min}$

- a spotting monitor feature and a CCD camera confirm actual spotting status on the PC screen
- supports micro-fractionation on up to 18 ABI MALDI-MS targets or 9 Bruker MALDI-MS targets (96 or 384 micro titer plate format)

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BioSpec-mini

The BioSpec-mini is a routine instrument with application-optimized user-interface designed particularly for the biosciences.

It is compact, easy-to-use and is also suitable for classic spectroscopic applications (190 - 1100 nm). It complements the UV-1240 compact spectrometer series and features a 5 nm fixed slit width.

- nucleic acid analysis, quantification of single and double stranded DNA, RNA

- Warburg-Christian method
- protein quantification according to, among others, the Lowry, Bradford and Biuret methods
- calculation of molecular mass, molar extinction coefficients as well as melting point determinations
- accessories for application-specific analyses (sippers, single and multi-cell holders, thermostatted cuvette holders, program cards [for instance kinetics], micro-cuvettes)

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
BIOTECH

To find your local Shimadzu contact
please visit www.shimadzu.eu

Shimadzu Europa GmbH
Albert-Hahn-Str. 6-10 · D-47269 Duisburg

Tel.: +49 - (0)203 - 76 87-0
Fax: +49 - (0)203 - 76 66 25
Email: shimadzu@shimadzu.de