

Preparative Supercritical Fluid Chromatography

Principles and applications

Shimadzu SFC User Meeting 2024

Isabelle François, PhD

Chromisa Scientific



Considerations in favor of SFC purification

All advantages of analytical SFC apply!

Advantages are even more significant at larger scale



Non-toxic
Pure, Stable



Non-flammable



Reclaimed from
Industrial plants



Easy fraction
recovery +
Limited solvent
evaporation



Profitable Total
Cost of
Ownership



Fast



Recyclable
Minimizes waste



Sustainable
Minimizes energy use




Applicable to many
small molecules




ATEX

No liquid
solvents
→ No ATEX





Sustainability and greenness across organisations



Our Corporate Commitments and operations to deliver

SYNGENTA GROUP


SEARCH

Innovation in

HOME / SUSTAINABILITY

We solve


Sustainability



Personal Care

Our Brands | Products and Applications | Technical Library | Discovery Zone | Sustainability | News | Contact

Progress beyond



Solvay > Sustainability

#Sustainable Development

Our sustainability strategy

Raising the bar to tackle climate change and resource scarcity to foster a better life

With Sustainable Development, Croda Corporation takes an integrated holistic approach, focusing on three pillars:

Sample basket | Login | English

Our Commitments

Animals

Non-communicable diseases

es in human and animal health

ty and safety

Boehringer Ingelheim

ities & Our People

ve and sustainable

Good Health for the Planet

Care for our environment



Sustainable chromatography tools

Green Chemistry Solvent Selection:

Preferred	Usable	Undesirable
Water	Cyclohexane	Pentane
Acetone	Heptane	Hexane(s)
Ethanol	Toluene	Di-isopropyl ether
2-Propanol	Methylcyclohexane	Diethyl ether
1-Propanol	Methyl t-butyl ether	Dichloromethane
Ethyl acetate	Isooctane	Dichloroethane
Isopropyl acetate	2-MethylTHF	Chloroform
Methanol	Tetrahydrofuran	Dimethyl formamide
Methyl ethyl ketone	Xylenes	N-Methylpyrrolidinone
1-Butanol	Dimethyl sulfoxide	Pyridine
T-Butanol	Acetic acid	Dimethyl acetate
	Ethylene glycol	Dioxane
		Dimethoxyethane
		Benzene
		Carbon tetrachloride

Greener Solvent Choices



- **SFC**
CO₂, Methanol, 2-Propanol
- **RP HPLC**
Acetonitrile, Water
- **NP HPLC**
DCM, Heptane, Ethyl Acetate





Tools to calculate “Greenness”

Analytical Method Greenness Score (AMGS)

Making the move towards modernized greener separations: introduction of the analytical method greenness score (AMGS) calculator†



Michael B. Hicks, ^{*a} William Farrell, ^b Christine Aurigemma, ^b Laurent Lehmann, ^c Lauren Weisel, ^a Kelly Nadeau, ^d Heewon Lee, ^e Carol Moraff, ^f Mengling Wong, ^g Yun Huang ^h and Paul Ferguson ⁱ



AMGS Calculator: <https://www.acsgcipr.org/amgs>

What is taken into account?

- Mass of solvent used and number of injections
- Safety, health and environmental impact of solvent
- Energy used during manufacture and disposal of solvent
- Energy used by instrument



Tools to calculate “Greenness”

Process Mass Intensity (PMI)

$$\text{Process Mass Intensity (PMI)} = \frac{\text{mass of all reactants}}{\text{mass of products}}$$




Each synthesis step increases the PMI with an average of 50 kg/kg



Greenness – SFC vs RPLC


Analytical Method Greenness Score (AMGS)

- ✓ *Lower solvent volumes*
- ✓ *Less waste, lower cost*



Equal
purification
success

64% less
waste
solvent per
sample



67% lower
running
cost per
sample

HPLC AMGS
1587


SFC AMGS
399

Data generated by Oncology Separation Sciences Team at Astra Zeneca
Calculated using Analytical Method Greenness Score (AMGS)



SFC purification in pharmaceutical industry

Pharmaceutical industry has been the main driver for SFC utilization, mainly due to the benefits in prep scale.






SFC scale – analytical, semi-prep, industrial/large



Waters UPC²
Analytical SFC



Waters SFC Prep100
Semi-preparative SFC




*Thar Process
Isolator SFC60*
**Large scale
preparative SFC**



SFC scale – analytical, semi-prep, industrial/large



Nexera Analytical SFC system
Analytical SFC
On-line SFE-SFC optional



Nexera Prep SFC system
Semi-preparative SFC



Providers of Prep SFC instrumentation

The providers listed below provide analytical and/or preparative packed-column SFC systems



smart with system



Teledyne LABS

ACCQ Prep
SFC



Waters

THE SCIENCE OF
WHAT'S POSSIBLE.™




Preparative Innovative Chromatography



(Very) large scale SFC systems are sometimes custom-built




Preparative SFC instrumentation





Preparative SFC instrumentation


What is different compared to LC purification?






Injection in Preparative SFC

Mixed stream injection



- Sample dissolved in the modifier, or other solvent
- Injection into the mixed mobile phase
- Significant impact of “strength” of injection solvent on peak shape at larger injection volumes
- Sample loop washed for a limited time
- Extensive wash of the injection port required

Modifier stream injection




- Principle used for preparative injections
- Sample dissolved in the modifier, or other solvent
- Dilution by CO₂ before the column without precipitation
- Applied at mobile phase strength
- Sample loop washed continuously
- Supporting overlapped injections in isocratic mode
- No sample carry over
- Only functions at sufficient modifier percentages (not below 5%)




Injection in Preparative SFC

Mixed stream injection



1. Sample dissolved in MeOH



Modifier stream injection



2. Sample dissolved in DMSO



Conditions: Sulfanilamide 5 mg/mL in MeOH or DMSO

Injection volume : 0.5, 1.0, 1.8 mL

Flow rate : 50 mL/min

Column : Diol 21.2 x 150 mm

BPR : 100 bar

Gradient : 20 to 50 % MeOH in CO₂ @ 10 %/min

T : 35°C




Injection in Preparative SFC


Stacked Injection Mode

- Increased throughput without compromising separation efficiency
- Only works in isocratic mode

Normal injections



Stacked injections







Injection in Preparative SFC

Stacked Injection Mode

- Example : hemp extract
- Analytical analysis was performed on the sample, to determine the appropriate conditions for purification
 - 5-40% gradient
- Method was scaled up
 - Isocratic method was developed @ 15% modifier
 - Three fractions




Courtesy of Chris Hudalla, ProVerde Laboratories, USA



Injection in Preparative SFC

Isocratic analytical method development for conditions to be transferred to SFC 200 for stacking experiments




SFC 200 Isocratic 30x150 mm

Condition Criteria

Cycle time: 3 min

Spacing: 1 min

Total elution time: 5 min






Fraction collection in Preparative SFC

- **CO₂ expansion after BPR makes collection more complicated**
 - Effluent is a biphasic mixture of CO₂ gas and micro-droplets (modifier and compound)
 - Equipment needed to separate gas from liquid without loss of material and without cross contamination
- **Heating is applied after the BPR to mitigate the cooling and associated issues (dry ice / blocked tubing)**
- **In many cases, a make-up solvent is added to avoid cross contamination and keep the purified compounds in solution (avoid blockage and precipitation). This only applies when the modifier percentage is low (below 5%)**






Fraction collection in Preparative SFC

Typically, cyclones or gas-liquid separators (GLS) are applied



Waters website

	Cyclone	GLS
Pressure & principle	High pressure Closed-bed format Liquid is pushed downwards due to design and gravimetry	Low pressure Open-bed format Spiral flow is created wherein the liquid and compounds are forced towards the wall due to centrifugal force. Spiralling downwards, a vortex is formed. When the diameter is sufficiently small, the clean CO ₂ moves upwards through the vortex and is sent to the top of the GLS.
Advantages 	No make-up required (CO ₂ remains its solvating power through the collection valve)	Cheaper materials can be applied Reduced safety risk Increased number of fractions
Disadvantages 	More expensive SS required due to high pressure Increased safety risk Limited number of fractions (dependent on the number of applied cyclones)	Additional make-up solvent required when using low modifier percentages Often fractions need to be combined Not suited for large quantity purification





Fraction collection in Preparative SFC

Shimadzu Nexera Prep SFC - LotusStream™



Shimadzu website

 **SHIMADZU**
Excellence in Science

New patented principle for gas-liquid separation :
LotusStream™

LotusStream separator (patented technology)

Decreases flow rate without increasing the pipe diameter by splitting flow through multiple channels. The CO₂ is discharged externally while the liquid travels along the column and drips directly into the collection vessel without dispersing or scattering the eluate.



**Chromisa
Scientific**



Fraction collection in Preparative SFC

Open bed vs Closed bed collection



Waters Prep SFC100

- **Open-bed**
- **Batch purification**
- **Smaller quantities / smaller campaigns / library purification**
- **Collection in tubes/vials/bottles**
- **1 GLS for all fractions**
 - **Sufficient solvent required (modifier / make-up) to avoid cross-contamination**
- **Detection : UV/PDA/ELSD/MS**




Thar Process Isolator SFC10

- **Closed bed**
- **Bulk purification**
- **Larger quantities / large campaign**
- **Collection in carboys**
- **Cyclone per fraction**
- **Detection : UV/ELSD**




MS & SFC orthogonality advantage

Mass-Directed Prep LC



Mass-Directed Prep SFC



Orthogonal
Compared to LC





Columns in Preparative SFC

- DAC : User packs column using a slurry of stationary phase mixed with solvent (MeOH, EtOH)
- Important considerations :
 - Availability and robustness of stationary phase in larger particles
 - Price of stationary phase
 - Eg DAC of 60 cm requires 65 kg stationary phase

Pre-packed



Hamilton website


Dynamic Axial Compression (DAC) column



Hanbon website



Columns in Preparative SFC - Robustness



These are hemp-derived products and there is inherent natural variability, but the chromatography shows repeatability and robustness.

Injection #19,684 on column
70mL of 50:50 distillate:ethanol
(10 months since column was packed)

CO₂ as mobile phase is more gentle than liquids for the stationary phase



Solvent recycling in Preparative SFC

- **For larger scale, CO₂ recycling is a must!**
 - **Improved sustainability and cost effectiveness**
- **Goal of recycler : liquidify “gaseous” CO₂ coming from collection system**



Preparative SFC – Screening and scaling


Requirements for high productivity purification

- High solubility in mobile phase
 - Large injection volumes → reduced resolution
 - Low solubility → poor peak shape (tailing, splitting)
 - Difficult to predict solubility in CO₂
- Short cycle time
- Preference to isocratic methods
- High stationary phase saturation capacity
 - Determines change in retention and peak shape with increasing sample load
- Good availability of stationary phase across particle size / prepacked columns / particles at affordable pricing
- High purity of fractions



Preparative SFC – Screening and scaling

- SFC uses compressible fluid as mobile phase
 - Particle size has significant impact on pressure drop
 - Pressure drop has significant impact on density and mobile phase strength
- When moving from analytical to prep, this pressure drop can be accounted for by adapting the BPR pressure
 - However, this quickly leads to unrealistic settings of BPR pressures and related overall pressures that the prep SFC systems can not achieve due system limitations
- Therefore, particle size is typically maintained from method development/screening to prep (5 or 10 μm)
- Fraction purity assessment can be done using smaller particles





Preparative SFC – Screening and scaling

- Following equations can be used while maintaining identical particle size in both analytical and prep SFC

$$F_{\text{Prep}} = F_{\text{Analytical}} * D_{\text{Prep}}^2 / D_{\text{Analytical}}^2$$

$$\text{Vol}_{\text{Prep}} = \text{Vol}_{\text{Analytical}} * D_{\text{Prep}}^2 / D_{\text{Analytical}}^2 * L_{\text{Prep}} / L_{\text{Analytical}}$$

Example:

Table 2: Strategy used to scale-up the separation when converting the CO₂ flow from mL/min to g/min.


	CO ₂ Flow	Co-solvent Flow
Method on UPC ²	1.275 mL/min	0.225 mL/min
CO ₂ mL to g conversion (density=0.936 g/mL)	1.1934 g/min	0.225 mL/min
Scaled to 19mm (19mm/3mm) ² x UPC ² Flow	47.87 g/min	9.025 mL/min
Final Prep Method	47.87 + 9.0205 = 56.9 g/min (Total Flow) 9.025/56.9*100=15.9% co-solvent	




Examples from the industry



Pharma API Purification
(Replacement for NPLC
and Orthogonal to RPLC)



Pharma enantiomer
separation
(chiral)



Synthetic Lipids
Purification



Natural Products
Extraction & Purification



Cannabis
Extraction &
Purification



Fish/algae Oil Omega
3-6 Fatty Acids
(DHA/EPA)



Examples from the industry

Fish/algae Omega 3 DHA/EPA

Multiple custom built large ID SFC systems (>50 cm ID)



KD-PUR* TECHNOLOGY

INNOVATION OF PURITY

● SUPERCritical FLUID TECHNOLOGY

Home > Supercritical fluid technology

SUPERCritical FLUID TECHNOLOGY

 **KD Pharma Group**[™]
Creating Health Solutions


<https://www.kdpharmagroup.com/en/our-difference/kd-pur-technology/supercritical-fluid-technology>



Examples from the industry

Cannabinoids

- THC remediation for hemp derived products
- THC/CBN/CBG isolation




Two Thar Process Isolator SFC60 systems (60 cm ID) + SFC10X as pilot scale system (10 cm ID) to assist method development



Beneficial Total-Cost-of-Ownership (TCO)

Example for cannabis purification and comparison to other techniques – solvent recycling and reduced evaporation of fractions and mobile phase






Examples from the industry

Lipid excipients

- Synthetic intermediate for vaccine drug research



Traditional purification method:

- Normal phase LC purification
- Toxic solvents are used – toluene, heptane, isopropyl acetate
- Large solvent consumption – up to 4,000 - 6,500L for 1 kg of material
- Consists of different steps of column conditioning (pre- and post run), gradient purification and column back-flushing – time consuming process
- Projects are under NDA, so not all details can be shared


CONFIDENTIAL



Examples from the industry

Lipid excipients

- Crude lipid starting at ~70% purity
- Purity criteria to meet specifications:
 - Total purity $\geq 97\%$
 - Individual impurities $< 0.5\%$









Examples from the industry

Lipid excipients

Approach R&D, Feasibility, scale-up



	<u>Phase I</u>	<u>Phase II</u>	<u>Phase III</u>	<u>Phase IV</u>
Goal	Analytical screening & <u>method development</u>	Semi-preparative purification	<u>Repeatability & loading study</u> for semi-preparative stage	<u>Scale-up to preparative fraction collection</u>
<u>Scale</u>	5-50 mg	mg – g		100 g - 10 kg
Column ID	4.6 mm	19 mm		10 cm
System	 <i>Thermo Vanquish incl CAD detection</i>  <i>Waters UPC²</i>	 <i>Waters SFC Prep 100 (mass directed)</i>		 <i>Thar Process Isolator™ SFC10X</i>



Examples from the industry

Lipid excipients

Column screening






Examples from the industry

Lipid excipients

Isocratic conditions on Isolator SFC10X




Courtesy of Thar Process



Examples from the industry

Lipid excipients – Ionizable lipid



Conditions:

Co-solvent: 7% methanol

Column: XXX (10 cm)

Cycle time: 21 minutes

~1.9 kg/day

Co-solvent usage : 230 L

Production of 1 kg

→ SFC up to 110 L

→ NPLC up to 6,500 L



Examples from the industry

Lipid excipients – Ionizable lipid – SFC vs NPLC

1 kg of purified product

SFC
110 L



NPLC
Up to
6,500 L

- Solvent recycling possible (80%)
→ 22 L consumed per 1 kg
- Only one solvent, no additives

- Mixture of solvents
- More challenging recycling and excessive solvent quantities



Acknowledgements:

- All contributors to this presentation
- AZ Mölndal for facilitation of this meeting
- Shimadzu & Gesa Schad
- You for your attention



Contact : Isabelle François
isabelle@chromisa.com
+32 478 999 726