


Oligonucleotide impurity analysis using an innovative SFC approach

Nov 7, 2024

Shinnosuke Horie, Shimadzu Europa GmbH


in place of Momoka Hayashida, Shimadzu Corporation



1. Introduction | Oligonucleotide therapeutics

Oligonucleotide therapeutics **draw** attention as one of the new drug modalities.

- Target at mRNA and miRNA that cannot be approached by existing drugs.
→ Expected to be **a treatment for diseases that have been difficult to treat.**
- 20 products were approved worldwide.
(except for mRNA therapeutics, as of June 2024).
(National Institute of Health Sciences, Division of Molecular Target and Gene Therapy Products, "Approved ON therapeutics," 2024-10-15)
- Composed of chemically modified oligonucleotides to improve biostability.
e.g. **phosphorothioate (PS) modification** was used in 12 approved oligonucleotide therapeutics.




PS modification

Quality control methods have been discussed


1. Introduction | Quality control issues

Co-elution of structurally related impurities



PS-modified oligonucleotide therapeutics and its impurities

A. Krotz *et al.*, *Bioorganic Med. Chem. Lett.* **1997**, 7, 73–78.



Chromatogram of synthesized oligonucleotide therapeutics (ion pair reversed-phased liquid chromatography; IPRP-LC)

J.M. Sutton. *et al.*, *J. Am. Soc. Mass Spectrom.* **2020**, 31(9), 1775-1782.


The method with different selectivity from IPRP-LC has been desired

➡ We focus on SFC

1. Introduction | Features of SFC with CO₂

Feature 1



High resolution



Nucleotide analysis: M.C. Beike *et al*, *J. Chromatogr. A* **2016**, 1436 84-90

Feature 2

Superior structural recognition




Peptide analysis: J. Molineau *et al.*, *J. Chromatogr. A* **2021**, 1658 462631

Feature 3

Low-polarity of CO₂

To the best of our knowledge, oligonucleotide analysis with SFC remains to be reported.

→controlling co-solvent compositions allows for the SFC applications to oligonucleotide.



Oligonucleotide
therapeutics

1. Introduction | Goal of this work

Evaluate the applicability of SFC separation to oligonucleotide analysis

①

Evaluation the separation performance for a 4-mer oligonucleotide, a model compound

- Optimization of column and modifier composition
- Separation performance for 4-mer with varying PS contents
- Retention behavior

(M. Hayashida *et al.*, *J of Chromatogr. A*, **2023**, 1708, 464333)

②

Apply to longer sequences


- Optimization of modifier composition for longer sequences
- Separation performance for deaminated products

(M. Hayashida *et al.*, *J of Chromatogr. A*, **2025**, 1744, 465731)

1. Materials | Analytes


Analytes		
Target	Sequence	[M-H] ⁻
a. T4	5'-T _o T _o T _o T-3'	1153.22
b. T4+1PS	5'-T _s T _o T _o T-3'	1169.20
c. T4+2PS	5'-T _s T _o T _s T-3'	1185.18
d. T4+3PS	5'-T _s T _s T _s T-3'	1201.15
e. TAGC	5'-T _o A _o G _o C-3'	1172.24
f. TAGC+1PS	5'-T _s A _o G _o C-3'	1188.22
g. TAGC+2PS	5'-T _s A _o G _s C-3'	1204.19
h. TAGC+3PS	5'-T _s A _s G _s C-3'	1220.17

T: thymidine, A: adenosine, G: guanosine, C: cytidine,
s :PS linkage, o:PO linkage



1-1. Results | Eluting condition

T4 were eluted from five types of column using the modifier containing ammonium formate



Analytical conditions

Flow rate : 1 mL/min

MP : CO₂/modifier = **20:80**

Make-up : 0.1 mL/min, MeOH

Column : Shim-pack UC
(4.6 X 150 mm; 3 μm)


BPR : 10.0 MPa, 50 °C

Oven temp. : 35 °C

Injection vol. : 1 μL, 100 μM

1-2. Results | 1st Column screening

We evaluate different types of column.



Analytical conditions

Flow rate : 1 mL/min

MP : CO₂/MeOH/1 M AmF
solution = **50:47.5:2.5**

Make-up : 0.1 mL/min, MeOH

Column : Shim-pack UC
(4.6 X 150 mm; 3 μ m)

BPR : 10.0 MPa, 50 °C


Oven temp. : 35 °C

Injection vol. : 1 μ L, 100 μ M



➔ **Diol II column** provides better separations.

1-3. Results | 2nd Column screening



We evaluate columns modified with polar stationary phases.




PVP




Triazole



NH₂



PBT



Name	m/z , [M-H] ⁻
■ T4	1153.22
■ T4+1PS	1169.20
■ T4+2PS	1185.18
■ T4+3PS	1201.15

*PS diastereomers

▲ Void volume (t_0)

Analytical conditions

Flow rate : 1 mL/min

MP : CO₂/MeOH/1 M AmF
solution = **50:47.5:2.5**

Make-up : 0.1 mL/min, MeOH

Column : Shim-pack UC
(**2.1** X 150 mm; 3 μ m)

BPR : 10.0 MPa, 50 °C



Oven temp. : 35 °C

Injection vol. : 1 μ L, 100 μ M

➔ **Diol II column** provides better separations.

1-4. Results | Optimize additives for hydrophilic sequences

- (e) TAGC
- (f) TAGC+1PS
- (g) TAGC+2PS
- (h) TAGC+3PS



Analytical conditions

Flow rate : **1.5** mL/min

MP : CO₂/additive in methanol buffer (95:5) = 50:50

Make-up : 0.1 mL/min, MeOH

Column : Shim-pack UC
(**2.1** X 150 mm; 3 μm)

BPR : 10.0 MPa, 50 °C


Oven temp. : 35 °C

Injection vol. : 1 μL, 100 μM


➔ **Aminoethanol acetate** provides better peak shape and separation performance

1-5. Results | Sequence-dependency for peak shape

- Good peak shapes were obtained in the sequences whose **GC content was less than 1/2**




N : theoretical plate numbers
 r : correlation coefficient




1-6. Results | Retention behavior

- Retention times of the selected sequences were **positively correlated with polar surface areas**



Hydrophilic interaction







Unique selectivity of Diol II
to PS contents

1. Summary

We evaluated the applicability of SFC to oligonucleotide analysis, using 4-mer oligonucleotides with different PS contents.

- Column screening revealed that Diol II provided better selectivity to PO and PS linkages.
- After examining several additives, we obtained good peak shapes and resolutions with aminoethanol.
- When the GC contents were less than 1/2, good peak shapes were obtained in the optimized method.



2. Introduction | Next topics

Evaluate the applicability of SFC separation to oligonucleotide analysis

①

Evaluation the separation performance
for a 4-mer oligonucleotide, a model compound

- Optimization of column and modifier composition
- Separation performance for 4-mer with varying PS contents
- Retention behavior

(M. Hayashida *et el.*, *J of Chromatogr. A*, **2023**, 1708, 464333)

②


Apply to longer sequences

- Optimization of modifier composition for longer sequences
- Separation performance for deaminated products

(M. Hayashida *et el.*, *J of Chromatogr. A*, **2025**, 1744, 465731)


2. Introduction | Deamination analysis

- The **4-amino group** of a cytidine, protected with acyl group, can be converted into a **hydroxyl group** due to the basic conditions of deprotection.
- Deamination analysis using MS is often challenging, and thus chromatographic separation is important.




A method with **different selectivity from ion pair reverse phase (IPRP) chromatography** has been desired.

2. Materials | Purification for chemically synthesized oligonucleotides



2-1. Results | Additives for the modifier



Analytical conditions


Flow rate: 1 mL/min
 Modifier: methanol/water=95:5
 containing 50 mM 2-aminoethanol or
 octylamine acetate
 B conc.: 20% (0–2 min), 20–60% (2–
 20 min), 60% (20–24 min), 60–20%
 (24–25 min), 20% (25–30 min)
 Make-up: 0.1 mL/min, methanol
 Column: Shim-pack UC-Diol II
 (4.6x150 mm; 3 μ m)
 BPR: 10.0 MPa, 50 $^{\circ}$ C
 Oven temp.: 35 $^{\circ}$ C
 Injection vol.: 1 μ L, 100 μ M

Peak	Sequence (DNA, 5'-3')	<i>m/z</i>
A	DMTr- TCACTTTTCAT	[M-2H] ²⁻ , 1626.32
A'	TCACTTTTCAT	[M-2H] ²⁻ , 1475.26
B	DMTr- TCACTTTTCATAATGCTGG	[M-4H] ⁴⁻ , 1440.51
B'	TCACTTTTCATAATGCTGG	[M-4H] ⁴⁻ , 1364.98

Peak	Sequence (5'-3', full 2'-MOE, C: 5-methyl cytosine)	<i>m/z</i>
1	DMTr- TCACTTTTCAT	[M-3H] ³⁻ , 1344.69
1'	TCACTTTTCAT	[M-3H] ³⁻ , 1243.97
2	DMTr- TCACTTTTCATAATGCTGG	[M-5H] ⁵⁻ , 1430.15
2'	TCACTTTTCATAATGCTGG	[M-5H] ⁵⁻ , 1369.73

- **Octylamine acetate** provided better peak shape for 10- and 18-mer
- Sharp peak for **2'-MOE-modified 18-mer** were detected

2-2. Results | Optimize additives



Peak	Sequence (5'-3', full 2'-MOE, <u>C</u> : 5-methyl cytosine)
1	DMTr-T <u>C</u> ACTTTTCAT
2	DMTr-T <u>T</u> ACTTTTCAT
3	DMTr-T <u>C</u> ACTTTTCATAATG <u>C</u> TGG
4	DMTr-T <u>T</u> ACTTTTCATAATG <u>C</u> TGG

Analytical conditions

Flow rate: 1 mL/min

Modifier: methanol/water=95:5 containing 50 mM additive and 50 mM acetic acid
 B conc.: 30% (0–5 min), 30–60% (5–30 min), 60% (30–35 min), 60–30% (36–40 min)

Make-up: 0.1 mL/min, methanol

Column: Shim-pack UC-Diol II (4.6x150 mm; 3 μm)


BPR: 10.0 MPa, 50 °C

Oven temp.: 35 °C

Injection vol.: 1 μL, 100 μM

- **Primary amine** provided better peak shapes than secondary and tertiary amines
- ➔ **Octylamine acetate** provided shaper peaks for 18-mer

2-3. Results | Optimize column oven temp.



Peak	Sequence (5'-3', full 2'-MOE, <u>C</u> : 5-methyl cytosine)
1	DMTr-T <u>C</u> ACTTTTCAT
2	DMTr-T <u>T</u> ACTTTTCAT
3	DMTr-T <u>C</u> ACTTTTCATAATG <u>C</u> TGG
4	DMTr-T <u>T</u> ACTTTTCATAATG <u>C</u> TGG

Analytical conditions

Flow rate : 1 mL/min

Modifier: methanol/water=95:5 containing 50 mM octylamine acetate

B conc.: 30% (0–5 min), 30–60% (5–30 min), 60% (30–35 min), 60–30% (36–40 min)

Make-up : 0.1 mL/min, methanol

Column : Shim-pack UC-Diol II (4.6x150 mm; 3 μm)

BPR : 10.0 MPa, 50 °C


Oven temp. : 35-50 °C

Injection vol. : 1 μL, 100 μM

50 and 60 °C provided better separation performance
→ 50 °C was selected

2-4. Results | Separate deaminated impurities contained in DMTr-on oligos

DMTr-on



Peak	Sequence(5'-3', full 2'-MOE, <u>C</u> : 5-methyl cytosine)
a	DMTr-TCAC T TTTCAT
b	DMTr-TT A CTTTTCAT
c	DMTr-TCAT T TTTCAT
d	DMTr-TCAC T TTT T AT
e	DMTr-TCAC T TTTCATAATGCTGG
f	DMTr-TT A CTTTTCATAATGCTGG
g	DMTr-TCAT T TTTCATAATGCTGG
h	DMTr-TCAC T TTT T AATAATGCTGG

Analytical conditions for SFC


Flow rate : 1.0 mL/min, MP : A)CO₂, B)50 mM octylamine acetate in methanol buffer (95:5) B conc.: 40% (10-mer), 45% (18-mer), Column : Shim-pack UC-Diol II (4.6 × 150 mm; 3 μm) BPR : 10.0 MPa, 50 °C, Oven temp. : 50 °C, Injection vol.: 1 μL, 100 μM

Analytical conditions for IPRP

Flow rate : 1.0 mL/min, MP : A)10 mM propylamine in water, B)ACN B conc.: 25% (10-mer), 23% (18-mer), Column : Shim-pack GIST C18 (4.6 × 150 mm; 3 μm), Oven temp. : 50 °C, Injection vol.: 1 μL, 100 μM

- **SFC** successfully separated some deaminated products from the 10- and 18-mer main products.
- **Sequence d** were coeluted with the main product using IPRP-LC.

2-5. Results | Separate deaminated impurities contained in DMTr-off oligos



Peak	Sequence(5'-3', full 2'-MOE, <u>C</u> : 5-methyl cytosine)
a'	TCACTTTCAT
b'	T T ACTTTCAT
c'	TCAT T TTTCAT
d'	TCACTTT T AT
e'	TCACTTTCATAATGCTGG
f'	T TACTTTCATAATGCTGG
g'	TCAT T TTTCATAATGCTGG
h'	TCACTTT T ATAATGCTGG
i'	TCACTTTCATAATG T TGG

Analytical conditions for SFC

Flow rate : 1.0 mL/min, MP : A)CO₂, B)50 mM octylamine acetate in methanol buffer (95:5) B conc.: 40% (10-mer), 45% (18-mer), Column : Shim-pack UC-Diol II (4.6 × 150 mm; 3 μm) BPR : 10.0 MPa, 50 °C, Oven temp. : 50 °C, Injection vol. : 1 μL, 100 μM

Analytical conditions for IPRP



Flow rate : 1.0 mL/min, MP : A)10 mM propylamine in water, B)ACN B conc.: 12% (10-mer), 13% (18-mer), Column : Shim-pack GIST C18 (4.6 × 150 mm; 3 μm), Oven temp. : 50 °C, Injection vol. : 1 μL, 100 μM

- DMTr-on and off sequences were almost **coeluted using SFC**.
- The separation capability **for DMTr-off sequences was higher** than that for DMTr-on sequences using IPRP-LC


2. Discussions | Different retention behavior

- The retention behavior of DMTr-on oligonucleotides using SFC was **dominated by interaction with oligonucleotides**.
- That using IPRP-LC was **strongly influenced by DMTr group** and thus the target sequences and their deaminated impurities were difficult to separate because of the similar hydrophobicity

➡ SFC could be a suitable method for analyzing oligonucleotides having **hydrophobic organ-targeting ligands** as well as the DMTr group



For more details...



2. Summary

- Octylamine acetate provided a sharp peak for 10- and 18-mer 2-MOE-modified oligonucleotides.
- 50 °C was better for column oven temp. to separate deaminated products.
- DMTr-on and –off oligos were co-eluted using SFC.
- Several deaminated impurities were separated from the main products using SFC.
- SFC and IPRP-LC exhibited **different** selectivity

1 & 2. Conclusion

This study demonstrated the applicability of SFC in oligonucleotide analysis

①

**Evaluation the separation performance
for a 4-mer oligonucleotide, a model compound**

- 4-mer oligonucleotides with different PS contents were separated

②

Apply to longer sequences

- We obtained peaks for 10- and 18-mer oligonucleotides
- Deaminated products were separated from the main products

Acknowledgement

- This research was a collaboration work with Prof. Satoshi Obika and Prof. Takao Yamaguchi (Graduate School of Pharmaceutical Sciences, Osaka University)
- This research was partially supported in part by AMED under Grant Number JP21ae0121022, JP21ae0121023, JP21ae0121024 (Project leader: Satoshi Obika).
- Co-authors:
Satoshi Obika and Takao Yamaguchi (Osaka University)
Junichi Masuda and Risa Suzuki (Shimadzu Corporation)
- Colleagues at Obika laboratory of Osaka University and Shimadzu Group
- Prof. Takeshi Bamba (Kyushu University)