

Nexera UC – Unified Chromatography

The latest addition to the chromatography toolbox

Dr. Gesa J. Schad
Shimadzu Europa GmbH




Nexera UC
Unified Chromatography

very A brief history of SFC

- Late 1800's: it was found that heavy, non-volatile organic compounds were surprisingly soluble in some inorganic gases above their critical point ('super'-critical).
- 1957: first ideas to use such gases as a chromatographic mobile phase
→ Acting as a solvent is a critical characteristic that differentiates the compressible mobile phase used in “**supercritical fluid chromatography**” (SFC) from gas chromatography (GC), where the mobile phase is considered to be an inert carrier
- Early 1990's: - concentration of polar modifiers was shown to be the primary retention control variable when added to the CO₂
 - pressure/density became a secondary control variable
 - use of additives allowed highly efficient elution of polar solutes, such as primary aliphatic amines → modern SFC was born

What is SFC and how does it work

- Supercritical CO₂ is a fluid state of carbon dioxide where it is held at or above its critical temperature (31.1 °C) and critical pressure (73.8 bar)




- Rule-of-thumb: any molecule that dissolves in methanol (or less polar solvents) is compatible with SFC
- CO₂ at its critical point is non-polar, solvent strength is increased by using a polar co-solvent
- When a co-solvent is introduced, the mobile phase is not truly supercritical, the terminology is used regardless

SFC vs. LC

- Equipment:

HPLC / UHPLC



SFC vs. HPLC

Similarities:

- RP-SFC, NP-SFC, IP-SFC and HILIC-SFC possible by changing polarity of stationary and mobile phase
- Organic modifiers and ionic additives can be used to adjust selectivity
- Packed column SFC is the most common methodology
- Silica based porous particles with various type of bonded phases are used for SFC separations
- Separations can be achieved isocratically or by using a solvent gradient

SFC vs. HPLC

Advantages of SFC:

- Low viscosity of the mobile phase
 - generates up to 5 times the linear velocity of RP-LC
 - less pressure drop across the column
- Resolving power is $\sim 5 \times$ that of HPLC
- Superior to LC for chiral separations
- Complimentary chromatographic selectivity to RP-LC
- “Green” technique – reduced organic solvent consumption
- High purity samples with little residual solvent obtained in preparative chromatography

SFE – what is it and how does it work


- spCO_2 has the transport capacity of a liquid, as the molecules are pressurized to give a fluid volume
- It has the mobility and penetrating power of a gas, as the molecules are capable of behaving independently
- spCO_2 can get into structures and surround molecules to be extracted using its gas-like properties and carry those analytes out of the matrix using its liquid-like properties.




solid



liquid



gas




supercritical fluid

Nexera UC – Features

- Fully automated online sample preparation and analysis
 - Faster, more efficient analytical work flow
 - Reduced quantitative error
 - Less degradation of labile samples
- Multi sample extraction by automated vessel changer
 - Max. 48 sample
- Flexible use as SFC or LC system
- Improved sensitivity due to low dead volume and low pulsation
BPR

Nexera UC – System Configuration


SFC / LC – UV System



- One Back Pressure Regulator after the detector
- No make-up pump needed
- Simple SFC / LC – UV hybrid system

Nexera UC – System Configuration

SFC / LC switching system – UV System




- One Back Pressure Regulator after the detector
- No make-up pump needed
- Simple SFC / LC – UV hybrid system

SFC or LC analysis in one batch



Easy batch file creation for LC/SFC switching using the Method Scouting Software




分析	バイアル番号	トレイ	サンプル名	メソッドファイル	データファイル	注入量
1	1	0		0151201_SFCLCtoSFC切り替え_tr.lcm	D:\20151201_SFCLCtoSFCtoLCData\2_004.lcd	1
2	-1	0		D:\20151201_SFCLCtoSFCtoLCBase.lcm	D:\20151201_SFCLCtoSFCtoLCData\2_005.lcd	1
3	1	0		D:\20151201_SFCLCtoSFCtoLCBase.lcm	D:\20151201_SFCLCtoSFCtoLCData\2_006.lcd	1
4	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ RP 0.1%AmAc, MeOH 5.40.012.lcd	1
5	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ RP 0.1%AmAc, MeOH 5.40.013.lcd	1
6	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ RP 0.1%AmAc, MeOH 5.40.014.lcd	1
7	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ RP 0.1%AmAc, MeOH 5.40.015.lcd	1
8	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ RP 0.1%AmAc, MeOH 5.40.016.lcd	1
9	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ RP 0.1%AmAc, MeOH 5.40.017.lcd	1
10	1	0		0151201_SFCLCtoSFCtoLC切り替え_tr.lcm	D:\20151201_SFCLCtoSFCtoLCData\2_014.lcd	1
11	-1	0		D:\20151201_SFCLCtoSFCtoLCBase.lcm	D:\20151201_SFCLCtoSFCtoLCData\2_015.lcd	1
12	1	0		D:\20151201_SFCLCtoSFCtoLCBase.lcm	D:\20151201_SFCLCtoSFCtoLCData\2_016.lcd	1
13	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ XR0DS, ACN, 0.1%FAw, 20.90.022.lcd	1
14	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ XR0DS, ACN, 0.1%FAw, 20.90.023.lcd	1
15	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ XR0DS, ACN, 0.1%FAw, 20.90.024.lcd	1
16	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ XR0DS, ACN, 0.1%FAw, 20.90.025.lcd	1
17	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ XR0DS, ACN, 0.1%FAw, 20.90.026.lcd	1
18	55	1	インダプロプライ	D:\20151201_SFCLCtoSFCtoLCBase.lcm	インダプロプライ XR0DS, ACN, 0.1%FAw, 20.90.027.lcd	1

SFC

LC

switching method LC → SFC is automatically inserted

SFC or LC analysis in one batch



LC


Mobile phase A : 0.1% Formic acid in H₂O

Mobile phase B : Acetonitrile

Flow rate: 1 mL/min

Gradient: 20 %B (0 min) – 90 %B (2-3 min)

Column: ShimPack XR-ODS II
75 x 3.0 mm, 2.2 µm



SFC

Mobile phase A: CO₂

Mobile phase B: 0.1 % Ammonium formate in MeOH


Flow rate: 1.5 mL/min

Gradient: 5%B(0min) – 90%B(2-3min)

Column: ShimPack UC-RP
150 x 2.1 mm, 3 µm

Nexera UC – System Configuration

Online SFC – PDA Scouting System



- One Back Pressure Regulator after the detector
- Two ovens with column switching valves for fast column screening for method development

Method scouting option

Method Scouting Solution for Nexera UC [Chiral]

New Open Save Save As Database System Check Info

Method Scouting Solution for Nexera UC

Estimated End Time : --- h --- min

Edit Monitor

A

☒ IA-3
☒ IB-3
☒ IC-3
☒ ID-3
☒ IE-3
☒ IF-3

Pump A

Pump B

☒ ACN
Approx. 0.1 L
☒ MeOH
Approx. 0.1 L
☒ ACN+EtOH
Approx. 0.1 L
☐
Approx. --- L

Total Inj Vol
0.04 mL

B

☒ AD-3
☒ AS-3
☒ AY-3
☒ OD-3
☒ OJ-3
☒ OZ-3

Method

Base Method File: C:\LabSolutions\Data\Project1\MSS_SFC_BaseMethod.lcm

Tray: 1 Inj. Vol.(uL): 1 # of Inj.: 1

Sample

Sample Name 1: Sample From 1 To 1
☐ Sample Name 2: From To
☐ Sample Name 3: From To
☐ Sample Name 4: From To
☐ Sample Name 5: From To

Data

Data Folder: C:\LabSolutions\Data\Project1

Filename: _ (Column)_ (PumpB)_ (Init.Conc.)_ (Final.Conc.)

Gradient

STOP: 12.00 min

Gradient Pattern Count 1

B.Conc(%)

40.0
20.0

6.00 0.00 2.00 4.00 min

Create Batch Add BatchQue

Nexera UC – System Configuration

Online SFE – SFC – MS System



- Additional extraction vessel
- Back Pressure Regulator to optimize SFE conditions
- Column switching valve for method development
- Make-up pump for ionization for MS detection

Extraction unit



SFE-30A

Extraction up to 80°C


Extraction vessels 0.2 mL or 5 mL



Large sample capacity



Extraction
Vessel




Up to **48 samples** can be automatically extracted and analysed by using the Rack Changer option

SFE sample preparation



Patented adsorbent effective for dehydration samples with high water content




Spot the sample



Cut DBS



Put DBS in extraction vessel



Extraction and analysis of 1 ppm phospholipid spiked into plasma on DBS

Sample preparation



Analysis of pesticides in spinach

Online SFE-SFC:

Extraction with CO₂ for 4 min (10 MPa, 40 °C)
(0 – 2 min static + 2 – 4 min dynamic extraction)

SFC Conditions:

Column: Inertsil ODS-EP HP, 4.6 × 250 mm
5 µm (GL Sciences Co.)

Modifier: 0.1% w/v (16 mM) NH₄OOCH in MeOH

Flow Rate: 3 mL/min

Gradient: 0 % (4 min); 11 % (15 min); 30 % (18min);
40 % (18.1 to 21 min)


BPR: 10 MPa

Detector: LCMS-8050 operating a MRM method

Make up: 0.1% w/v (16 mM) NH₄OOCH in MeOH


Flow Rate: 0.1 mL/min

Spinach (spiked 100ppb pesticides)






Nexera UC – Features

- High sensitivity detection
 - Improved sensitivity due to low dead volume BPR
 - Splitless injection into the MS without band broadening




x 5




Nexera UC – Features

- High sensitivity detection
 - Concentrated sample due to evaporation of CO₂ after BPR



x 6



LC Conditions

Column : FC ODS, 2.1×150 mm, 3μm
 Mobile Phase A : 0.1% W/V Ammonium Formate in MeOH
 Mobile Phase B : 0.1% W/V Ammonium Formate in H₂O
 B Conc. : 30%
 Flow Rate : 0.4 mL/min
 Detector : LCMS-8040
 sample : prostaglandin 100 μg/L



SFC Conditions

Column : Inertsil Diol, 2.1 × 150 mm, 3 m
 Modifier : 0.1 %w/v ammonium formate in MeOH
 Conc. : 20 %
 Flow Rate : 1 mL/min
 BPR : 10 MPa (50 °C)
 Detector : LCMS-8040
 sample : prostaglandin 100 μg/L

Application

Analysis of coenzyme Q₁₀ in supplements

→ Conventional Method: Coenzyme Q₁₀ is prone to oxidation




→ No oxidation with SFE

Extraction with 5 % MeOH in CO₂ for 4 min

Gradient: 5 %B (0 min – 4 min for SFE)

5 % - 50 %B in 5 min






Application

Extraction of polymer additives

→ Conventional Method: 10 – 20 h soxhlet extraction

→ SFE: 7 min CO₂ extraction

- Crushed polymer samples are filled in the extraction vessel
- After 7 min CO₂ extraction the sample is ready for SFC analysis



Application

Cleaning validation



→ Conventional Method:

- Sampling – solvent extraction – concentration – HPLC analysis



→ SFE-SFC:

- Swab is enclosed in the extraction vessel for SFE





Application

Chiral Method Scouting

→ Dedicated Method Scouting Solution:

- Selection of up to 12 columns
- Solvent blending of modifier



Column: CHIRALPAK AD-H
Modifier: Acetonitrile/ Methanol=7/3
Resolution = 2.66 (Rt = 2.644)

Thank you for your attention !