

Determination of heavy metals in flavored e-liquids using ICP-OES Spectrometry

Ludivine Fromentoux<sup>1</sup>, Uwe Oppermann<sup>2</sup>, Jan Knoop<sup>2</sup>, and Mirna Markusi<sup>3</sup>  
<sup>1</sup>Shimadzu France SAS, Marne La Vallee, France, <sup>2</sup>Shimadzu Europa GmbH, Spectroscopy, Duisburg, Germany  
<sup>3</sup>Shimadzu d.o.o. Croatia, Zagreb, Croatia

Overview

Between 2013 and 2015, the number of electronic cigarettes users increased by 86%. Excluding USA and China, seven countries including France, Germany, Poland, Italy and Russia; represent more than 75% of the e-cigarettes world market.



Figure 1: electronic cigarette: safe pleasure or heavy metals poison ?

In order to prevent European consumers health, the European Council included some regulations for this kind of product in the Tobacco Products Directive (TPD) 2014/40/EU. Concerning heavy metals contamination in e-liquids, European Commission directive 88/388/CEE is applied. This regulation sets maximum concentration values for 5 toxic elements : arsenic, cadmium, mercury, lead and tin.

Element	As	Cd	Hg	Pb	Sb
Maximum authorized Value (mg/L)	3	1	1	10	5

Table 1: Heavy metals maximum authorized values in e-liquids according directive 88/388/CEE.

A new French experimental standard XP D90-300-2 was recently developed, forcing the e-liquids producer to quantify this five elements in their products. Moreover, recent studies show that aluminum, copper, iron, chromium and nickel could also be found in e-liquids and e-vapor.

Then a fast and easy to use method for simultaneous analysis of 10 elements (Al, As, Cd, Cr, Cu, Fe, Hg, Ni, Pb, Sb) in e-liquids has been developed using the Shimadzu ICPE-9820. This method is also able to quantify the 10 elements in e-liquids raw materials (polypropylene-glycol PG and Vegetable Glycerin VG) and e-vapors.

Shimadzu ICPE-9820: performance, simplicity and reduced costs.



The combination of a large CCD detector (1 inch) and an efficient vacuum optical system in ICPE-9820 results in high sensitivity and resolution (5 pm at 200 nm). The Shimadzu exclusive vertical Minitorch limits memory effects and drastically reduce Argon consumption during analysis run (12 L/min). Moreover, the new Shimadzu RF solid state generator insure a great stability for plasma generation and measurement. At the same time, this device enables an “eco mode” which reduces Argon total consumption at less than 4 L/min between measurement run. For each sample, the easy to use ICPEsolution software saves all the wavelength between 167 nm and 800 nm, thereby facilitating the data processing.

Samples preparation and method parameters

Four different flavored e-liquids, called ELiQ1 to ELiQ4 are analyzed. One of them, ELiQ2, is studied with two different nicotine level (3 and 18 mg/L). Tow raw materials, vegetable glycerine (VG) and propylene glycol (PG) are also studied. Finally, e-vapors from ELiQ4 are analyzed, thanks to a smoking machine (50 and 180 puffs).

Thanks to ICPE-9820 robustness and performance, only easy sample dilution is used. Only 1.5 mL of nitric acid 15% are added at 2 g of sample. Then, the resulting solution is completed at 25 mL with distilled water before ICP analysis.

Parameter	Setting
RF generator power	1.2 kW
Plasma gas	10 L/min
Auxilliary gas	0.6 L/min
Carrier gas	0.7 L/min
Nebulizer type	Concentric 10
Spray chamber	Cyclonic
Plasma torch	Minitorch
View	axiale
Exposure time	30 s

Table 2: ICPE-9820 relevant parameters

Calibration

For each studied element, calibration curves including 4 or 6 points in the targeted concentration range are realized in 2% nitric acid solutions. All samples are measured in triplicate and one of them, ELiQ1 is spiked with known concentration values (50 ppb) in order to check the method accuracy. A 1 ppm Yttrium internal standard solution in 2% nitric acid was mixed using online mixing with sample, before it is aspirated.




Figure 2: Calibration curves obtained using 2% HNO<sub>3</sub>. Red elements are specifically targeted by 88/388/CEE.

The different curves in Figure 2, show all correlation coefficients r<sup>2</sup> achieving a 0.9999 level.

Detection limits calculation

Detection limits (LD) are calculated automatically by ICPEsolution software with 3σ method. All LD shown in table 3 below are reported in e-liquids sample.

Element	As	Cd	Hg	Pb	Sb
LD (ppb)	150	1	75	33	54
Element	Al	Cr	Cu	Fe	Ni
LD (ppb)	7.5	2.5	2.5	6.3	6.3

Table 3: LD values for e-liquids sample.

As shown in Table 3, all the LD obtained values are far lower than the maximum authorized contents in e-liquids required by the European regulation. It demonstrates ICPE-9820 ability to check the safety of this kind of product.

Results

The results of five studied e-liquids are collected in Table 4. Table 5 refers to element values in raw materials and e-vapors.

	ELiQ1	ELiQ2 (3mg)	ELiQ2 (18 mg)	ELiQ3	ELiQ4
Al (ppm)	0.062	0.091	0.097	0.051	< LD
As (ppm)	0.025	0.020	0.012	0.080	0.069
Cd (ppm)	< LD	< LD	< LD	3.7.10 <sup>-3</sup>	7.2.10 <sup>-3</sup>
Cr (ppm)	< LD	< LD	< LD	3.0.10 <sup>-3</sup>	2.6.10 <sup>-3</sup>
Cu (ppm)	7.5.10 <sup>-3</sup>	0.046	0.049	< LD	0.020
Fe (ppm)	0.105	0.094	0.035	< LD	0.366
Hg (ppm)	< LD	< LD	< LD	< LD	< LD
Ni (ppm)	< LD	< LD	< LD	< LD	5.3.10 <sup>-3</sup>
Pb (ppm)	0.076	0.244	0.256	0.275	0.265
Sb (ppm)	0.109	0.363	0.394	0.346	0.410

Table 4: Elements concentration in studied e-liquids.

The quantitation results in table 4 demonstrate that ICPE-9820 is able to quantify simultaneously the various target elements present in e-liquid samples whatever the flavor or the nicotine level.

	PG	VG	ELiQ4 vapor 50 puffs	ELiQ4 vapor 180 puffs
Al (ppm)	0.097	0.086	0.028	0.043
As (ppm)	<LD	0.055	0.02	<LD
Cd (ppm)	0.012	0.017	<LD	<LD
Cr (ppm)	0.012	0.014	1.6.10 <sup>-3</sup>	1.7.10 <sup>-3</sup>
Cu (ppm)	<LD	<LD	0.020	0.011
Fe (ppm)	0.151	0.014	0.035	0.078
Hg (ppm)	<LD	<LD	<LD	<LD
Ni (ppm)	7.7.10 <sup>-3</sup>	0.016	3.0.10 <sup>-3</sup>	3.5.10 <sup>-3</sup>
Pb (ppm)	0.268	0.219	0.01	0.049
Sb (ppm)	0.384	0.324	0.028	0.038


Table 5: Elements concentration in studied e-liquids raw materials and ELiQ4 vapours.

As demonstrated by the values in Table 5, the ICPE-9820 using the developed method is also able to quantify element concentrations in e-liquids raw materials and vapors.

Method accuracy

In order to demonstrate the method accuracy, some spike of each element is done in e-liquid ELiQ1. The achieved results are mentioned in graph 1 and calculated according to:

recovery(%) = (value after spike - initial value) / initial value x 100



Graph 1: recovery values for each element in ELiQ1.

The recovery rate for all the elements are between 80% and 120% (min. 87%, max. 115%) indicating the suitability of the developed method from point of accuracy, independent from initial element concentration.