

Analysis of antiparasitic agents using LC-MS

Veterinary pharmaceuticals including antibiotics and hormones are used to prevent disease in livestock, promote growth, and enhance the feed efficiency. Antiparasitic agents are also widely used to eliminate parasites from the alimentary canal. Residual standards are being established for antiparasitic agents as residual levels in meat present similar health problems to antibiotics and hormones.

The four components used for this test were 5-hydroxythiabendazole, thiabendazole, flubendazole, and albendazole. Their structures are shown in Fig. 1. A residual standard is set for each of these components in food. HPLC is prescribed for the

analysis of these components but LC-MS permits analysis with extremely high selectivity and sensitivity.

Fig. 2 shows the LC-MS analysis results for the four antiparasitic agents. Each component could be positively identified using mass chromatography at the mass number of the protonated molecule of each component. Fig. 3 shows their mass spectra. Selected ion monitoring (SIM) permits highly sensitive analysis.

Fig. 4 shows the calibration curves in the range from 10 to 1000 ppb. Each curve shows good linearity.




Fig. 1 Structures of antiparasitic agents






Fig. 2 UV absorption (305nm) and mass chromatograms of antiparasitic agents


5-hydroxythiabendazole



Thiabendazole



Flubendazole



Albendazole




Fig. 3 Mass spectra of antiparasitic agents

Fig. 4 Calibration curves(10ppb-1000ppb)

Table 1 Analytical conditions for LC-MS

| | |
|---------------------|--|
| Column | : Inertsil ODS-2 (2.1 mmI.D. x 150 mm) |
| Mobile phase A | : 5mM acetic acid ammonium acetate buffer |
| Mobile phase B | : acetonitrile |
| Gradient program | : 0%B (0 min) - 100%B (20 min) |
| Flow rate | : 0.2 mL/min |
| Column temperature | : 40 °C |
| Probe voltage | : +4.5 kV (ESI-Positive mode) |
| Nebulizing gas flow | : 4.5 L/min |
| CDL voltage | : -40 V (0 - 14 min), -50 V (14.01 - 25 min) |
| DEFs voltage | : +40 V (0 - 14 min), +45 V (14.01 - 25 min) |
| Scan range | : m/z 100-400 |

SHIMADZU CORPORATION International Marketing Division

3. Kanda-Nishikicho 1-chome, Chiyoda-ku, Tokyo 101-8448, Japan

Phone: 81 (3) 3219-5641 Fax: 81 (3) 3219-5710

Cable Add. SHIMADZU TOKYO