

Rapid Screening and Quantitation of Pesticide Residues in *Cannabis* by Modified QuEChERS and LC-MS-MS

ASMS 2016 TP 230

Jeffrey H. Dahl¹, Julie Kowalski², Derek Laine³,
and Jason Zitzer³

¹Shimadzu Scientific Instruments; Columbia, Maryland

²Restek Corporation; Bellefonte, Pennsylvania

³Trace Analytics; Spokane, Washington

Rapid Screening and Quantitation of Pesticide Residues in *Cannabis* by Modified QuEChERS and LC-MS-MS

Introduction

Medical and recreational use of marijuana (*Cannabis* spp.) is expanding in the United States at a rapid pace, and domestic production has increased more than ten-fold in the last 25 years. This extremely high value crop is vulnerable to mold and insects so growers frequently apply pesticides and antifungals to protect their plants. These chemical residues may pose a danger to

consumers, so highly sensitive and selective methods for their detection in the complex cannabis matrix are required. We developed a rapid and effective LC-MS-MS method with modified QuEChERS sample preparation for detection of nearly 200 chemical residues in dried cannabis flower and used the method to test a wide selection of products offered for commercial sale.




Figure 1 Modified QuEChERS Extraction

Method

Test portions of dried cannabis flower were homogenized and extracted using a modified QuEChERS extraction with dispersive SPE cleanup (Restek cat no. 26237 and 26243). Grinding was performed with a SPEX GenoGrinder. Detection was carried out by LC-MS-MS using a Shimadzu Prominence HPLC with LCMS-8050 triple

quadrupole mass spectrometer. Electrospray ionization was used with continuous polarity switching to measure analytes in both modes throughout the run. Pesticide recovery was determined using spiking experiments and matrix-matched calibration curves. All analysis was carried out in a Washington state certified cannabis testing lab.

Table 1 LCMS-8050 Instrument parameters

LC Column	: Restek ARC-18 (2.1×100 mm, 3 μm)
Mobile Phase A	: 5 mM NH4OAC + 0.1% Formic Acid
Mobile Phase B	: Methanol
LC Flow Rate	: 0.5 mL/min
Heater Gas	: 10 L/min
Interface Temp	: 400 °C
Nebulizing Gas	: 3 L/min
Drying Gas	: 10 L/min
DL Temp	: 250 °C
Heat Block Temp	: 400 °C

Rapid Screening and Quantitation of Pesticide Residues in *Cannabis* by Modified QuEChERS and LC-MS-MS



Figure 2 Cannabis samples for pesticide testing




Figure 3 Representative chromatogram and number of compounds by recovery (inset)

Results and Discussion

QuEChERS extraction with dispersive SPE cleanup provided the best combination of pesticide recovery and cleanup for dried flower cannabis samples. Matrix matched calibration curves were linear within the quantitation limits established for each compound, which was compound dependent, but ranged from as low as 20 ppb or lower to greater than 500 ppb for a few substances. Detection limits and quantitation limits were required to have 3:1 and 10:1 signal to noise respectively, and quantitation limits were required to have less than 20% RSD in triplicate analyses. Recovery was compound

dependent however the majority were within the range of 70-120% while outliers above and below the range were observed. The method was validated in three different cannabis strains using matrix matched calibration curves and triplicate QC spikes at three levels. In a subset of randomly tested cannabis flower samples offered for commercial sale, the three most commonly detected residues were piperonyl butoxide, myclobutanil, and boscalid. Concentrations for pesticides ranged from the detection limit up to the microgram/gram level.

Rapid Screening and Quantitation of Pesticide Residues in *Cannabis* by Modified QuEChERS and LC-MS-MS




Figure 4 Representative individual chromatograms at the 100 ng/g dried cannabis spiking level (Abamectin chromatogram from 1 mcg/g level)

Rapid Screening and Quantitation of Pesticide Residues in *Cannabis* by Modified QuEChERS and LC-MS-MS

Table 2 Detection of pesticides in 39 dried flower samples offered for retail sale

Survey of pesticides in dried <i>cannabis</i> flower offered for retail sale					
Sample	Residue detected	mcg/g	Sample	Residue detected	mcg/g
A			P	Piperonyl butoxide	0.32
B			Q	Imidacloprid	0.49
C			R		
D			S	Piperonyl butoxide	0.69
E			T		
F			U	Myclobutanil	0.02
G	Piperonyl butoxide	0.53	V		
H	Piperonyl butoxide	0.05	W	Piperonyl butoxide	12.46
I			X	Piperonyl butoxide	0.16
J	Spinosyn A	0.02		Fipronil	0.04
	Spinosyn D	0.02	Y		
K	Myclobutanil	0.01	Z		
L	Myclobutanil	0.02	AA		
M	Dinotefuran	13.44	AB		
	Boscalid	81.79	AC	Piperonyl butoxide	14.99
	Pyraclostrobin	0.40	AD	Permethrin	0.35
	Trifloxystrobin	0.08	AE		
	Fludioxonil	0.42	AF	Piperonyl butoxide	0.08
	Myclobutanil	1.21	AG		
N	Dinotefuran	1.20	AH		
	Boscalid	5.79	AI		
	Fludioxonil	0.01	AJ	Diuron	0.06
	Myclobutanil	0.08	AK	Piperonyl butoxide	0.13
O	Boscalid	0.15	AL		
	Myclobutanil	0.05	AM		

Rapid Screening and Quantitation of Pesticide Residues in *Cannabis* by Modified QuEChERS and LC-MS-MS




Figure 5 Representative calibration curves

Rapid Screening and Quantitation of Pesticide Residues in *Cannabis* by Modified QuEChERS and LC-MS-MS



Figure 6 Sample preparation

Table 3 Pesticides detected in 39 dried flower samples

Number of pesticide detections in dried <i>cannabis</i> flower		
Residue name	Detections	Rate
Boscalid	3	7.7%
Dinotefuran	2	5.1%
Diuron	1	2.6%
Fipronil	1	2.6%
Fludioxonil	2	5.1%
Imidacloprid	1	2.6%
Myclobutanil	6	15.4%
Permethrin	1	2.6%
Piperonyl butoxide	9	23.1%
Pyraclostrobin	1	2.6%
Spinosyn A	1	2.6%
Spinosyn D	1	2.6%
Trifloxystrobin	1	2.6%
One or more	19	49%

Rapid Screening and Quantitation of Pesticide Residues in *Cannabis* by Modified QuEChERS and LC-MS-MS




Figure 7 RSD for each compound in triplicate QC samples at the 50 ng/g spike level, in three matrices

Conclusion

A validated method for detection of chemical residues in dried cannabis flower samples was developed. Our method can detect low levels of common pesticides in samples offered for retail sale with excellent selectivity and speed. Measurements of a larger selection of commercially available cannabis samples are being carried out.

First Edition: June, 2016



Shimadzu Corporation
www.shimadzu.com/an/

For Research Use Only. Not for use in diagnostic procedure.

This publication may contain references to products that are not available in your country. Please contact us to check the availability of these products in your country.

The content of this publication shall not be reproduced, altered or sold for any commercial purpose without the written approval of Shimadzu. Company names, product/service names and logos used in this publication are trademarks and trade names of Shimadzu Corporation or its affiliates, whether or not they are used with trademark symbol "TM" or "@". Third-party trademarks and trade names may be used in this publication to refer to either the entities or their products/services. Shimadzu disclaims any proprietary interest in trademarks and trade names other than its own.

The information contained herein is provided to you "as is" without warranty of any kind including without limitation warranties as to its accuracy or completeness. Shimadzu does not assume any responsibility or liability for any damage, whether direct or indirect, relating to the use of this publication. This publication is based upon the information available to Shimadzu on or before the date of publication, and subject to change without notice.