

# Application News

Liquid Chromatograph Mass Spectrometer LCMS-8060

## Analysis of Antidepressants in Serum / Plasma Using RECIPE® ClinMass® TDM Kit System with Fully Automated Sample Preparation LC/MS/MS System

Anja Grüning<sup>1</sup>, Akihiro Kunisawa<sup>2</sup>, Ionela Regos<sup>3</sup>

1 Shimadzu Europa GmbH, 2 Shimadzu Corporation, 3 RECIPE Chemicals + Instruments GmbH

### User Benefits

- ◆ Full solution provided by Shimadzu and RECIPE®
- ◆ Fully automated sample preparation
- ◆ Verified method for RECIPE® ClinMass® TDM Kit System for Antidepressants in Serum / Plasma

### Introduction

Major depressive disorder (MDD) affects many adults, often treated initially in primary care with medications. Standard treatments include classical tricyclic antidepressants (TCAs) and newer options like selective serotonin reuptake inhibitors (SSRIs), serotonin-norepinephrine reuptake inhibitors (SNRIs), and tetracyclic antidepressants (TeCAs). Therapeutic drug monitoring (TDM) is essential for optimizing doses and minimizing side effects, particularly in cases of uncertain adherence, suboptimal tolerability, or drug interactions.

RECIPE's fully validated analytical method provides the quantification of 37 Antidepressants and metabolites (Table 2) for TDM using LC-MS/MS.<sup>[1]</sup>

By addition of the Shimadzu CLAM (Clinical Laboratory Automated sample preparation Module) in front of the LC-MS/MS system (Figure 1) the required sample preparation could be fully automated which achieves results on a fast and high-precision analytical workflow.

To prove that the automated sample preparation leads to reliable results a method verification procedure was evaluated according to the CLSI Guidelines EP06-A, EP15-A3, EP17-A2.

Then the samples were loaded directly into the CLAM-2040. It was programmed to perform protein precipitation using Precipitant P including internal standards from the ClinMass® TDM Kit System for Antidepressants followed by filtration and sample collection. The sample is then transported using an arm from the CLAM-2040 to the LC without human intervention for LC-MS/MS analysis. Due to overlapped sample preparation (Figure 2) and analysis the throughput was one complete analysis each 5.0 min. Analytical conditions are listed in Table 1. The optimized MRM transitions are summarized in Table 2.




Fig. 2 Scheme fully automated sample preparation and analysis



Fig. 1 CLAM LCMS TQ

### Materials and Methods

Fast, sensitive and robust LC-MS/MS systems provide the basis for routine analysis in clinical laboratories. For the described verification, a Shimadzu CLAM-2040 coupled with a Nexera X3 UHPLC system and a LCMS-8060 triple-quadrupole mass spectrometer was used.

Thirty-seven Antidepressants in serum were verified using the ClinMass® TDM Platform (order no. MS9000) in combination with the ClinMass® Add-on Set for Antidepressants in Serum / Plasma (order no. MS9400) (RECIPE®, Germany).

Lyophilized, matrix-based calibrator and control samples were reconstituted, aliquoted and stored until use.

Table 1 Analytical conditions

Mass Spectrometer	: LCMS-8060
Ionization	: Electrospray Ionization (ESI), positive
Interface Voltage	: 4 kV
Heating Gas	: 10 L/min
DL Temp.	: 250 °C
Interface Temp.	: 300 °C
Nebulizing Gas	: 3 L/min
Drying Gas	: 10 L/min
Heat Block	: 400 °C
CID	: 270 kPa
UHPLC	: Nexera X3
Column Oven	: 40 °C
Injection Volume	: 2.0 µL
Flow rate	: 0.7 mL/min
Time Programme	: Binary gradient

Time [min]	Mobile Phase A [%]	Mobile Phase B [%]
Initial	95	5
0.10	95	5
0.20	75	25
1.50	50	50
2.50	45	55
2.60	20	80
3.00	20	80
3.10	95	5
3.70	95	5

Table 2 MRM transitions and parameters of the analytes and isotope-labelled substances

Analyte / IS	Quantifier MRM		Dwell Time msec	Q1 Pre Bias [V]	CE [V]	Q3 Pre Bias [V]
	Precursor [m/z]	Product [m/z]				
<b>Agomelatine</b>	244	185.1	5	-27	-27	-19
<b>Atomoxetine</b>	256.1	117.1	5	-24	-35	-11
<b>Bupropion</b>	240	184.1	5	-16	-12	-12
<b>Citalopram</b>	325.2	109.1	5	-12	-15	-11
<b>Clomethiazole</b>	162	126.2	5	-28	-20	-24
<b>Desmethylcitalopram</b>	311	109.1	5	-11	-10	-11
<b>Desmethylfluoxetine</b>	296	134.1	5	-26	-7	-22
<b>Desmethylnianserine</b>	251.2	208.1	5	-10	-20	-21
<b>Desmethylnortazapine</b>	252.2	195.1	5	-17	-35	-19
<b>Desmethylsertraline</b>	291.9	158.9	5	-29	-25	-29
<b>Dosulepin</b>	296	202	5	-14	-40	-21
<b>Duloxetine</b>	298	154	5	-19	-6	-15
<b>Dihydro-Bupropion*</b>	242	168.1	5	-17	-3	-30
<b>Fluoxetine</b>	310.1	148.1	5	-17	-9	-15
<b>Fluvoxamine</b>	319	71.2	5	-16	-33	-13
<b>Guanfacine</b>	245.9	60	5	-12	-16	-23
<b>Hydroxybupropion</b>	256	139	5	-17	-40	-14
<b>Methylphenidate</b>	234	84	5	-23	-21	-16
<b>Mianserine</b>	265	58	5	-30	-23	-22
<b>Milnacipran</b>	247.2	100.1	5	-29	-30	-10
<b>Mirtazapine</b>	266	195.1	5	-19	-15	-20
<b>Moclobemide</b>	269	182.1	5	-18	-40	-19
<b>Nefazodone</b>	470.1	274.2	5	-17	-40	-28
<b>O-Desmethyltramadol</b>	250.1	58.1	5	-40	-5	-40
<b>O-Desmethylvenlafaxine</b>	264.2	107	5	-30	-32	-10
<b>Opi Pramol</b>	364.1	171.1	5	-11	-40	-17
<b>Paroxetine</b>	330	70.1	5	-10	-29	-12
<b>Reboxetine</b>	314	176.1	5	-22	-5	-25
<b>Ritalinic acid</b>	220	84.1	5	-23	-40	-14
<b>Sertraline</b>	306.1	159	5	-15	-15	-25
<b>Tianeptine</b>	437	292	5	-16	-35	-20
<b>Tramadol</b>	264.2	58	5	-30	-5	-21
<b>Tranlycypromine</b>	134	117.1	60	-10	-14	-23
<b>Trazodone</b>	372	120.1	5	-24	-54	-21
<b>Venlafaxine</b>	278.1	58.1	5	-13	-10	-23
<b>Vortioxetine</b>	299	150	5	-20	-24	-15
<b>d3-Agomelatine</b>	247	185.1	5	-17	-16	-12
<b>d3-Atomoxetine</b>	259.2	117	5	-10	-13	-21
<b>d9-Bupropion</b>	249.1	185.1	5	-10	-12	-21
<b>d6-Citalopram</b>	331.2	109.1	5	-35	-40	-22
<b>d3-Dosulepin</b>	299.2	203.1	5	-11	-35	-18
<b>d7-Duloxetine</b>	305.3	154.1	5	-15	-6	-18
<b>d3-Desmethylcitalopram</b>	314.2	109	5	-11	-22	-22
<b>d5-Desmethylfluoxetine</b>	301.1	139.1	5	-11	-7	-29
<b>d4-Desmethylsertraline</b>	295.9	159.9	5	-11	-28	-18
<b>d5-Fluoxetine</b>	315.1	153.1	5	-12	-8	-29
<b>d6-O-Desmethyltramadol</b>	256.2	64	5	-40	-5	-40
<b>d3-Fluvoxamine</b>	322.2	74.1	5	-12	-32	-14
<b>d6-Hydroxybupropion</b>	262.1	139	5	-10	-35	-16
<b>d9-Methylphenidate</b>	243.2	93	5	-27	-21	-16
<b>d3-Mianserine</b>	268.2	61.2	5	-10	-24	-24
<b>d10-Milnacipran</b>	257.2	110.1	5	-18	-30	-20
<b>d3-Mirtazapine</b>	269.2	195.1	5	-10	-15	-22
<b>d8-Moclobemide</b>	277	183.1	5	-10	-30	-21
<b>d6-Nefazodone</b>	476.2	280.1	5	-15	-29	-20
<b>d4-Opi Pramol</b>	368.1	175.1	5	-25	-16	-19
<b>d4-Paroxetine</b>	334.2	74.1	5	-10	-30	-14
<b>d5-Reboxetine</b>	319.2	176.1	5	-12	-5	-13
<b>d3-Sertraline</b>	309.1	159	5	-15	-25	-18
<b>d6-Tramadol</b>	270.1	64	5	-30	-5	-21
<b>d9-threo-Dihydro-Bupropion</b>	251.1	169	5	-16	-25	-19
<b>d5-Tranlycypromine</b>	139.1	122.1	60	-16	-34	-25
<b>d6-Venlafaxine</b>	284.2	58.2	5	-29	-34	-24
<b>d8-Vortioxetine</b>	307.1	153	5	-11	-24	-17

Table 2 MRM transitions and parameters of the analytes and isotope-labelled substances (continued)

Analyte / IS	Qualifier MRM		Dwell Time msec	Q1 Pre Bias [V]	CE [V]	Q3 Pre Bias [V]
	Precursor [m/z]	Product [m/z]				
<b>Agomelatine</b>	244	141.1	5	-27	-47	-14
<b>Atomoxetine</b>	256.1	65	5	-24	-55	-13
<b>Bupropion</b>	240	131.1	5	-16	-25	-13
<b>Citalopram</b>	325.2	262.1	5	-12	-20	-18
<b>Clomethiazole</b>	162	70.9	5	-28	-44	-17
<b>Desmethylcitalopram</b>	311	234.1	5	-11	-26	-24
<b>Desmethylfluoxetine</b>	296	105.2	5	-26	-23	-19
<b>Desmethylmianserine</b>	251.2	118.1	5	-10	-25	-23
<b>Desmethylmirtazapine</b>	252.2	209.1	5	-17	-23	-22
<b>Desmethylsertraline</b>	291.9	275	5	-29	-12	-13
<b>Dosulepin</b>	296	178.1	5	-14	-47	-18
<b>Duloxetine</b>	298	97.1	5	-19	-53	-19
<b>Dihydro-Bupropion*</b>	242	116.1	5	-17	-20	-11
<b>Fluoxetine</b>	310.1	65.1	5	-17	-80	-14
<b>Fluvoxamine</b>	319	200	5	-16	-20	-22
<b>Guanfacine</b>	245.9	123.2	5	-12	-46	-21
<b>Hydroxybupropion</b>	256	238.1	5	-17	-3	-40
<b>Methylphenidate</b>	234	56	5	-23	-50	-20
<b>Mianserine</b>	265	77.1	5	-30	-55	-14
<b>Milnacipran</b>	247.2	230.2	5	-29	-22	-15
<b>Mirtazapine</b>	266	72.1	5	-19	-21	-13
<b>Moclobemide</b>	269	139	5	-18	-10	-13
<b>Nefazodone</b>	470.1	246.1	5	-17	-35	-16
<b>O-Desmethyltramadol</b>	250.1	232.1	5	-40	-10	-40
<b>O-Desmethylvenlafaxine</b>	264.2	58	5	-30	-5	-22
<b>Opipramol</b>	364.1	143.1	5	-11	-45	-14
<b>Paroxetine</b>	330	192.1	5	-10	-21	-20
<b>Reboxetine</b>	314	91	5	-22	-20	-19
<b>Ritalinic acid</b>	220	56.1	5	-23	-46	-20
<b>Sertraline</b>	306.1	275.1	5	-15	-18	-18
<b>Tianeptine</b>	437	228	5	-16	-38	-15
<b>Tramadol</b>	264.2	42	5	-30	-55	-15
<b>Tranlycypromine</b>	134	65	60	-16	-34	-25
<b>Trazodone</b>	372	78.1	5	-24	-30	-14
<b>Venlafaxine</b>	278.1	260.2	5	-13	-12	-18
<b>Vortioxetine</b>	299	109	5	-20	-40	-22

\*Dihydro-bupropion is determined as the sum of erythro- and threo-dihydro-bupropion

## ■ Results

The trueness was determined by 4-fold analysis of two different quality control (QC) samples in a single analysis sequence. The results (precision in CV% and deviation from the target in % Bias) are summarized in Table 3.

To determine the precision two different levels of QC samples were prepared in 8-fold and analyzed in a single analysis sequence. The intraassay precision for each level is summarized in Table 4.

For determination of the linearity and the lower limit of quantification (LLOQ) several dilutions of the ClinCal® Serum Calibrator Set lyophil. for Antidepressants (order no. MS9413, RECIPE®, Germany) were prepared in 3-fold and analyzed in a single analysis sequence. The results for linearity evaluation and for the LLOQ are summarized in Table 5.

Table 3 Trueness of measurement

Analyte	Target value	Measured value	CV	Bias	Target value	Measured value	CV	Bias
	[µg/L] MS9482 Level I	[µg/L] Mean (n=4)			[µg/L] MS9482 Level II	[µg/L] Mean (n=4)		
<b>Agomelatine</b>	32.9	29.3	3.2	-10.9	315	313	1.8	-0.6
<b>Atomoxetine</b>	392	405	9.0	3.4	917	953	5.2	4.0
<b>Bupropion</b>	39.4	41.6	3.5	5.5	92.5	98.4	1.6	6.4
<b>Citalopram</b>	44.6	43.3	2.9	-2.9	104	108	2.9	3.5
<b>Clomethiazole</b>	450	482	4.7	7.1	2764	2840	2.4	2.8
<b>Desmethylcitalopram</b>	51.6	49.4	2.1	-4.2	121	114	1.5	-6.0
<b>Desmethylfluoxetine</b>	124	114	3.6	-7.9	292	282	2.7	-3.6
<b>Desmethylmianserine</b>	28.3	31.6	2.9	11.5	69.1	76.9	2.5	11.4
<b>Desmethylmirtazapine</b>	37.8	41.0	2.6	8.5	86.0	92.0	3.5	6.9
<b>Desmethylsertraline</b>	39.2	39.4	2.0	0.4	87.3	94.9	5.5	8.8
<b>Dosulepin</b>	44.3	42.0	3.3	-5.2	103	97.5	2.0	-5.3
<b>Duloxetine</b>	53.0	54.9	2.0	3.7	123	127	2.1	3.5
<b>Dihydro-Bupropion*</b>	268	260	1.4	-3.2	627	622	3.0	-0.8
<b>Fluoxetine</b>	101	99.4	3.7	-1.6	237	239	3.5	1.0
<b>Fluvoxamine</b>	109	101	2.8	-7.4	242	245	2.2	1.4
<b>Guanfacine</b>	2.62	2.39	6.3	-8.9	6.26	6.0	1.1	-4.1
<b>Hydroxybupropion</b>	418	414	1.9	-0.9	972	977	2.0	0.5
<b>Methylphenidate</b>	9.26	9.20	2.4	-0.6	21.7	21.7	3.1	0.0
<b>Mianserine</b>	32.4	31.2	2.0	-3.9	76.3	74.2	0.8	-2.7
<b>Milnacipran</b>	65.1	62.6	3.8	-3.8	153	160	5.2	4.7
<b>Mirtazapine</b>	35.5	34.4	0.7	-3.2	83.4	79.8	2.1	-4.3
<b>Moclobemide</b>	497	464	3.4	-6.6	1041	1059	2.8	1.7
<b>Nefazodone</b>	105	99.5	1.0	-5.2	250	231	1.8	-7.4
<b>O-Desmethyltramadol</b>	233	206	2.2	-11.6	501	487	4.2	-2.7
<b>O-Desmethylvenlafaxine</b>	109	107	2.7	-2.0	251	252	1.6	0.3
<b>Opipramol</b>	109	105	2.6	-3.8	258	253	1.6	-2.1
<b>Paroxetine</b>	42.1	40.8	2.8	-3.0	97.9	99.9	1.9	2.0
<b>Reboxetine</b>	143	134	2.3	-6.5	330	320	2.9	-3.1
<b>Ritalinic acid</b>	70.7	72.0	2.7	1.9	161	168	4.4	4.3
<b>Sertraline</b>	29.2	30.0	1.9	2.7	152	152	2.1	-0.2
<b>Tianeptine</b>	29.0	31.0	1.3	6.7	70.8	73.8	2.4	4.3
<b>Tramadol</b>	237	229	1.7	-3.2	529	521	3.4	-1.5
<b>Tranlycypromine</b>	19.4	18.9	7.8	-2.7	47.8	46.6	3.7	-2.5
<b>Trazodone</b>	555	547	2.5	-1.4	1247	1305	1.8	4.7
<b>Venlafaxine</b>	61.6	61.7	1.5	0.1	142	149	2.1	5.1
<b>Vortioxetine</b>	22.8	22.6	0.8	-0.8	49.8	52.6	1.2	5.7

\*Dihydro-bupropion as a sum of erythro and threo-dihydro-bupropion

Table 4 Intraassay results [CV%]

Analyte	Sample	Measured value	CV	Sample	Measured value	CV
		[µg/L]	[%]		[µg/L]	[%]
		Mean (n=8)	[%]		Mean (n=8)	[%]
<b>Agomelatine</b>	MS9482 Level I	30.3	4.4	MS9482 Level II	311	1.8
<b>Atomoxetine</b>	MS9482 Level I	395	7.4	MS9482 Level II	947	3.6
<b>Bupropion</b>	MS9482 Level I	41.2	2.6	MS9482 Level II	98.5	2.1
<b>Citalopram</b>	MS9482 Level I	43.7	2.5	MS9482 Level II	108	2.0
<b>Clomethiazole</b>	MS9482 Level I	475	4.9	MS9482 Level II	2823	1.9
<b>Desmethylcitalopram</b>	MS9482 Level I	48.3	3.7	MS9482 Level II	116	2.9
<b>Desmethylfluoxetine</b>	MS9482 Level I	115	3.4	MS9482 Level II	283	3.3
<b>Desmethylmianserine</b>	MS9482 Level I	31.3	2.8	MS9482 Level II	75.4	3.2
<b>Desmethylmirtazapine</b>	MS9482 Level I	41.1	2.3	MS9482 Level II	91.6	2.8
<b>Desmethylsertraline</b>	MS9482 Level I	39.0	3.1	MS9482 Level II	93.7	4.5
<b>Dosulepin</b>	MS9482 Level I	42.1	3.0	MS9482 Level II	98.5	1.9
<b>Duloxetine</b>	MS9482 Level I	54.0	4.1	MS9482 Level II	127	2.1
<b>Dihydro-Bupropion*</b>	MS9482 Level I	260	1.4	MS9482 Level II	617	2.2
<b>Fluoxetine</b>	MS9482 Level I	101	3.5	MS9482 Level II	240	3.1
<b>Fluvoxamine</b>	MS9482 Level I	99.9	2.2	MS9482 Level II	243	2.4
<b>Guanfacine</b>	MS9482 Level I	2.47	6.5	MS9482 Level II	6.08	2.6
<b>Hydroxybupropion</b>	MS9482 Level I	410	2.6	MS9482 Level II	974	1.4
<b>Methylphenidate</b>	MS9482 Level I	9.21	2.3	MS9482 Level II	21.6	2.7
<b>Mianserine</b>	MS9482 Level I	31.0	1.4	MS9482 Level II	74.3	0.6
<b>Milnacipran</b>	MS9482 Level I	63.6	3.7	MS9482 Level II	159	4.4
<b>Mirtazapine</b>	MS9482 Level I	34.1	2.2	MS9482 Level II	80.5	2.0
<b>Moclobemide</b>	MS9482 Level I	467	2.6	MS9482 Level II	1055	2.8
<b>Nefazodone</b>	MS9482 Level I	98.9	1.1	MS9482 Level II	231	1.4
<b>O-Desmethyltramadol</b>	MS9482 Level I	206	2.3	MS9482 Level II	492	4.1
<b>O-Desmethylvenlafaxine</b>	MS9482 Level I	106	2.9	MS9482 Level II	249	1.6
<b>Opipramol</b>	MS9482 Level I	105	2.2	MS9482 Level II	250	2.5
<b>Paroxetine</b>	MS9482 Level I	41.2	3.4	MS9482 Level II	98.5	2.1
<b>Reboxetine</b>	MS9482 Level I	132	2.5	MS9482 Level II	318	2.1
<b>Ritalinic acid</b>	MS9482 Level I	73.0	3.2	MS9482 Level II	167.2	3.9
<b>Sertraline</b>	MS9482 Level I	29.7	2.2	MS9482 Level II	150.8	2.3
<b>Tianeptine</b>	MS9482 Level I	31.2	2.1	MS9482 Level II	72.7	2.5
<b>Tramadol</b>	MS9482 Level I	231	2.7	MS9482 Level II	525	2.7
<b>Tranlycypromine</b>	MS9482 Level I	19.3	7.6	MS9482 Level II	47.0	4.0
<b>Trazodone</b>	MS9482 Level I	550	2.1	MS9482 Level II	1322	1.8
<b>Venlafaxine</b>	MS9482 Level I	61.5	1.1	MS9482 Level II	149	1.7
<b>Vortioxetine</b>	MS9482 Level I	22.3	1.7	MS9482 Level II	52.6	1.6

\*Dihydro-bupropion as a sum of erythro and threo-dihydro-bupropion

Table 5 Linearity evaluation, including LLOQ / LOD

	Linear Range [µg/L]	R <sup>2</sup>	LLOQ [µg/L]	LOD [µg/L]
Agomelatine	4.95 - 683	0.999	4.95	1.65
Atomoxetine	69.0 - 2065	0.992	69.0	23.0
Bupropion	3.40 - 224	1.000	3.40	1.13
Citalopram	6.75 - 233	1.000	6.75	2.25
Clomethiazole	110 - 6336	0.995	110	36.7
Desmethylcitalopram	8.45 - 271	0.997	8.45	2.82
Desmethylfluoxetine	20.1 - 672	0.996	20.1	6.7
Desmethylmianserine	10.7 - 185	0.998	10.7	3.57
Desmethylmirtazapine	12.4 - 194	0.993	12.4	4.13
Desmethylsertraline	12.6 - 199	0.999	12.6	4.2
Dosulepin	3.80 - 234	0.999	3.80	1.27
Duloxetine	17.5 - 274	0.991	17.5	5.83
Dihydro-Bupropion*	25.5 - 1578	0.999	25.5	8.50
Fluoxetine	31.2 - 535	0.995	31.2	10.4
Fluvoxamine	31.7 - 527	0.999	31.7	10.6
Guanfacine	0.834 - 13.7	0.998	0.834	0.278
Hydroxybupropion	130 - 2075	0.996	130	43.3
Methylphenidate	0.705 - 47.8	0.999	0.705	0.235
Mianserine	5.20 - 173	1.000	5.20	1.73
Milnacipran	22.4 - 397	0.993	22.4	7.47
Mirtazapine	5.65 - 179	0.999	5.65	1.88
Moclobemide	53.3 - 2691	0.997	53.3	17.8
Nefazodone	7.40 - 465	0.999	7.40	2.47
O-Desmethyltramadol	73.9 - 1185	0.994	73.9	24.6
O-Desmethylvenlafaxine	18.9 - 579	0.999	18.9	6.3
Opi Pramol	8.80 - 590	0.997	8.80	2.93
Paroxetine	12.6 - 220	0.997	12.6	4.2
Reboxetine	23.6 - 752	0.999	23.6	7.87
Ritalinic acid	13.6 - 382	0.991	13.6	4.53
Sertraline	4.70 - 293	0.999	4.70	1.57
Tianeptine	2.88 - 165	0.992	2.88	0.96
Tramadol	75.4 - 1138	0.998	75.4	25.1
Tranlycypromine	6.15 - 108	0.999	6.15	2.05
Trazodone	144 - 2552	0.994	144	48.0
Venlafaxine	20.6 - 349	1.000	20.6	6.87
Vortioxetine	0.751 - 118	0.999	0.751	0.250

\*Dihydro-bupropion as a sum of erythro and threo-dihydro-bupropion

## ■ Conclusion

The ClinMass® TDM Kit System for Antidepressants in Serum / Plasma (order no. MS9000 and MS9400) was successfully verified on the CLAM-2040 with the analytical system LCMS-8060 from Shimadzu.

All 37 analytes passed the acceptance criteria for accuracy (trueness, precision) and linearity.

The lower limit of quantification (LLOQ) was below published clinical reference ranges.

## ■ References

1. Instruction Manual, ClinMass® TDM Kit System, Antidepressants in Serum / Plasma, RECIPE® Chemicals + Instruments GmbH



Shimadzu Corporation  
www.shimadzu.com/an/

SHIMADZU Europa GmbH,  
www.shimadzu.eu