



REACH  
**SEPARATIONS**

## **Sustainability considerations for small scale purification workflows - SFC or HPLC?**

---

Stéphane Dubant


Shimadzu SFC User Meeting  
October 2025

# Reach Separations

- Established over 2 sites
- Experts in **Purification**
- Provide several **analytical** services
- Recently acquired by the **Catsci** Group




Nottingham




Strasbourg



# What we do at Reach



**Historically  
for purification**





**The Reach  
Paradigm**



# Green Focus – *Energy*

## Energy consumption of each purification stage:

- Analytical scale
- Preparative scale
- Dry-Down



| AMGS score |     |
|------------|-----|
| HPLC       | SFC |
| 59         | 58  |
| 1575       | 411 |

➡ The **drying** process = most **energy intensive** stage

➡ LC utilises more than **twice the energy** as SFC

*Based on manufacturer Data & actual measures in the lab*

**1 kWh = 87 g eq. CO<sub>2</sub>**

Average value in France in 2022  
(Source: *electricity maps*)

# Green Focus - AGMS

| AGMS score        | NP    | SFC  |
|-------------------|-------|------|
| Analytical scale  | 22.8  | 10.8 |
| Preparative scale | 504.4 | 85.9 |



**Biggest** impact on **preparative** scale



**SFC** significantly improves **green score**



# Chiral Prep example

| AMGS calculator comparison |                            |         |               |                |           |          |             |           |            |
|----------------------------|----------------------------|---------|---------------|----------------|-----------|----------|-------------|-----------|------------|
| Nb of injections           | Nb of analytes of interest | Diluent | Column        | Mobile phase A | Flowrate  | Run-time | % cosolvent | Cosolvent | AMGS score |
| 10                         | 2                          | MeOH    | 30x250mm, 5µm | CO2            | 150mL/min | 5 min    | 20          | MeOH      | 330.56     |
| 10                         | 2                          | MeOH    | 30x250mm, 5µm | CO2            | 150mL/min | 5 min    | 20          | EtOH      | 735.27     |
| 10                         | 2                          | MeOH    | 30x250mm, 5µm | CO2            | 150mL/min | 5 min    | 20          | iPOH      | 512.68     |
| 10                         | 2                          | MeOH    | 30x250mm, 5µm | CO2            | 150mL/min | 5 min    | 20          | MeCN      | 4110.3     |
| 7                          | 2                          | MeOH    | 30x250mm, 5µm | Heptane        | 42mL/min  | 5min     | 20          | EtOH      | 910.14     |

|    |   |      |               |     |           |       |   |    |        |
|----|---|------|---------------|-----|-----------|-------|---|----|--------|
| 10 | 2 | MeOH | 30x250mm, 5µm | CO2 | 150mL/min | 5 min | 0 | NA | 122.91 |
|----|---|------|---------------|-----|-----------|-------|---|----|--------|

Solvent cost in our lab (1L or 1kg):  $\text{CO}_2 < \text{MeOH} < \text{iPOH} < \text{EtOH} < \text{MeCN}$

$\text{CO}_2$  -> 7tons tank outside the building (food grade)

Solvents -> come in 30L shuttle drum (prep HPLC grade)



# Why is MeOH considered greener than EtOH

## 1. Lower Cumulative Energy Demand

Methanol generally requires less energy to produce and purify compared to ethanol, especially when considering synthetic routes from biomass or CO<sub>2</sub>. This lower energy footprint contributes positively to the AMGS score. [\[About the... – ACSGCIPR\]](#)

## 2. Solvent Health and Safety Profile

While both solvents are flammable and toxic to some extent, methanol is often used in smaller volumes in chromatography and has a well-understood risk profile. Ethanol, although less toxic, may have higher exposure risks due to its volatility and broader use in larger volumes.

## 3. Environmental Impact

Methanol biodegrades relatively quickly and has a lower potential for bioaccumulation. It also produces fewer harmful byproducts during combustion or disposal. Ethanol, while biodegradable, is often derived from agricultural sources, which can raise concerns about land use and sustainability.

[\[Methanol v...inability?\]](#)

## 4. Waste and Instrument Efficiency

Methanol is compatible with many high-efficiency chromatographic methods (like UHPLC and SFC), which use less solvent and generate less waste.


Ethanol, due to its viscosity and polarity, may require longer run times or higher volumes, increasing waste and energy usage. [\[Analytical...Calculator\]](#)

## 5. Feedstock Versatility

Methanol can be produced from a wide range of feedstocks including biomass, CO<sub>2</sub>, and even waste gases, making it more adaptable to green production methods. Ethanol is primarily derived from crops, which can compete with food production and require significant water and land resources



# Pressure based method development in SFC



Pressure gradient profile

In this case, 100 to 350 bar (instrument limitation).

Then method optimisation, typically if we can, we go isobaric.

We are working without co-solvent:

- Only the stationary phase will affect selectivity (large set needed)
- Eluting strength tweaking is done by **adjusting CO<sub>2</sub> density**



# Adjusting CO<sub>2</sub> density

There are 2 ways you can affect the density of the fluid in the system:

- Change the pressure in the system:
  - Change the BPR value (100 to 200 bar)
  - Change the flow-rate (2 to 3mL/min)
- Change the temperature

If Pressure



then density



, typically compounds elute faster

If Temperature




then density



, typically compounds elute later



# Playing on BPR pressure




Fish oil extract

As pressure increases, compounds are eluting faster.



# Playing on flow-rate



Harder to spot as decreasing flow-rate obviously changes retention time but it will reduce pressure on the head of the column therefore compounds are typically more retained.

When Flowrate




, density



typically compounds elute **later!**

Mix: Amino-Biphenyl, Benzophenone, Dibromobiphenyl

# Playing on temperature



When Temperature



, density




, typically compounds elute **later!**




# Transfer from analytical to analytical

System with lots of valves  
and narrow tubings




Simple system, larger  
tubings



<sup>1</sup> Mix: Amino-Biphenyl, Benzophenone, Dibromobiphenyl



# Transfer from analytical to analytical



Elution order changed!

Column Head  
Pressure: 180Bar

$\Delta P=80$  bar

BPR Pressure:  
100Bar


Column Head  
Pressure: 120Bar

$\Delta P=20$  bar


BPR Pressure:  
100Bar

The BPR pressure value is important but so is pressure drop!

# Fragrance compounds



# Method « repeatability »



Natural product extract

Overlay of 5 injections in pressure gradient mode




# Semi volatile target in plant extract



# Semi volatile target from a plant extract

Product is an oil and contains about 24% of the taret molecule (by GC-FID)  
This is for indsutrial scale SFC (tons/years)




Current process is based on many successive distillations to go from 24% to 95% purity (GC-FID)

Signal: FID1A

| RT [min] | Type | Area   | Area% | Name   |
|----------|------|--------|-------|--------|
| 10.078   | BB   | 1482.1 | 52.3  |        |
| 10.864   | VB   | 380.5  | 13.4  |        |
| 11.853   | VV   | 36.5   | 1.3   |        |
| 12.750   | VB   | 693.8  | 24.5  | Target |



# Semi volatile target from a plant extract



Full screening of columns/conditions  
Optimised method: CO<sub>2</sub>/EtOH 97/3 on a BiP column.

The method was then scaled to 3cm ID column.  
Oil injected neat.



# Semi volatile target from a plant extract


## Input

| Crude (g) | Target content (GC-FID) | Maximum recoverable Target (g) |
|-----------|-------------------------|--------------------------------|
| 15,5      | 24,50%                  | 3,7                            |

## Output

| Recovered Target (g) | Purity (GC-FID) | Recovery |
|----------------------|-----------------|----------|
| 3,4                  | 91,20%          | 81%      |

Successful project, recovery is good (after dry-down) considering this is a semi-volatile compound! (BP around 110°C)



| Signal: FID1A |        |       |        |
|---------------|--------|-------|--------|
| RT [min]      | Area   | Area% | Name   |
| 9.978         | 25.2   | 2.0   |        |
| 10.788        | 8.1    | 0.6   |        |
| 11.769        | 13.4   | 1.1   |        |
| 12.075        | 7.7    | 0.6   |        |
| 12.557        | 1.7    | 0.1   |        |
| 12.600        | 5.8    | 0.5   |        |
| 12.685        | 1162.9 | 91.2  | Target |



# Sustainability considerations




|                          |         |                  |         |
|--------------------------|---------|------------------|---------|
| Method Number:           | Percent | Greenness Score: | Percent |
| 2025-09-06-22:11:02.109  |         | 441.23           |         |
| Instrument Energy Score: | 293.76  | 66.58%           |         |
| Solvent Energy Score:    | 9.95    | 2.25%            |         |
| Solvent EHS Score:       | 137.52  | 31.17%           |         |

More difficult to compare but this method allows the removal of about 9 distillations steps (customer feedback), so likely it is an overall improvement in terms of sustainability





# RP analytical data



Target is 1.6% in UV at 480nm  
Identity was confirmed by RT, MS and UV spectrum.  
50-98% MeCN gradient



# SFC analytical screening approach

A set of 10 columns were tested to find the best selectivity ( all in 4.6x250mm, 10µm).

A pressure gradient approach was tested but did not give enough eluting strength (need for a ***higher-pressure system!***)

A solvent gradient (ethanol -> customer requirement) approach has then been used.


Best column from the screening set was Venusil HILIC.

The provided sample was dissolved in DCM (not good but the only solvent that could dissolve everything...).


Analysis was performed with ELSD and UV (480nm) and prep was performed at 480nm.



# Promising screening data




The ELSD allows us to see that the non-UV content of the sample elutes before the UV-visible one.




This combination of column and solvent seems promising.



# Optimised Isocratic SFC method





# Prep chromatogram example



Retention time was very stable over the course of 40 injections.



# Profile crude/fraction comparison



| # | Name        | RT (min) | Area (mAU-s) | Area%  | Height (mAU) | Height% | UV Conf. Match | FacSymmetry | Width (min) |
|---|-------------|----------|--------------|--------|--------------|---------|----------------|-------------|-------------|
| 1 |             | 1.548    | 3.045        | 0.18   | 1.023        | 0.15    |                | 0.80973     | 0.098       |
| 2 | Astaxanthin | 3.328    | 1382.317     | 81.9   | 600.056      | 86.35   | 1000           | 0.81788     | 0.242       |
| 3 |             | 3.557    | 279.825      | 16.579 | 87.877       | 12.65   |                | 0.72995     | 0.169       |
| 4 |             | 5.163    | 8.578        | 0.508  | 2.781        | 0.4     |                | 0.87974     | 0.103       |
| 5 |             | 5.537    | 6.656        | 0.394  | 1.936        | 0.28    |                | 0.94291     | 0.109       |
| 6 |             | 6.1      | 7.382        | 0.437  | 1.272        | 0.18    |                | 0.98843     | 0.203       |



# Outcome

- We have been able to find SFC conditions that seem suitable for astaxanthin purification from the given matrix
- The mobile phase contains only CO<sub>2</sub> & ethanol
- The UV purity (@480nm, absorption maximum of astaxanthin and its family) goes from 1.6 to 80% with 1 pass on SFC
- Actual published methods are:
  - Low pressure NP (with DCM & Acetone)
  - High Pressure LC (with high content of ACN)
  - CPC (*n*-hexane–ethanol–water)


|                          |         |                  |  |
|--------------------------|---------|------------------|--|
| Method Number:           |         | Greenness Score: |  |
| 2025-09-07-21:43:38.556  |         | 1749.16          |  |
| Instrument Energy Score: | 1175.04 | 67.18%           |  |
| Solvent Energy Score:    | 38.53   | 2.20%            |  |
| Solvent EHS Score:       | 535.59  | 30.62%           |  |



# Small Synthetic peptide




# Customer sample





| Retention Time (min) | Area   | Area% |
|----------------------|--------|-------|
| 2.25                 | 14.6   | 0.41  |
| 2.28                 | 1655.9 | 46.72 |
| 2.40                 | 22.4   | 0.63  |
| 2.69                 | 78.5   | 2.21  |
| 2.95                 | 1773.3 | 50.03 |

RP analysis with formic acid  
PEP5 is composed of N amino acids to separate from a N-1 amino-acids one.

Peak RT: 2.30



# SFC conditions screening



## *Columns tested:*


- Venusil HILIC
- Polar RP
- Gemini C8
- Prep PhenHex
- Luna NH2

4.6x250mm, 10µm

Several mobile phase tested  
MeOH/H<sub>2</sub>O 95/5 Amm For 20mM



# SFC loading study



3mg of crude diluted into 600 $\mu$ L of MeOH (good solubility)


Injection volumes: 5, 10, 20 & 40 $\mu$ L

Resolution still looked ok at 40 $\mu$ L so selected for the prep

Shimadzu SFC User Meeting - October 2025




# Preparative run example




11 injections of 40 $\mu$ L were performed so 2.2mg on column total (200 $\mu$ g/injection)

# Results



| Retention Time (min) | Area   | Area%  |
|----------------------|--------|--------|
| 2.96                 | 1033.9 | 100.00 |

|        | Amount (mg) | UV Purity (%) |
|--------|-------------|---------------|
| PEP5_1 | 0.9         | 98            |
| PEP5_2 | 0.8         | 84            |



| Retention Time (min) | Area   | Area% |
|----------------------|--------|-------|
| 1.97                 | 21.3   | 1.78  |
| 2.26                 | 17.3   | 1.44  |
| 2.29                 | 1010.3 | 84.31 |
| 2.70                 | 108.5  | 9.06  |
| 2.98                 | 40.9   | 3.42  |

Elution order is reversed between RP & SFC



# AMGS considerations

|                                                          |        |
|----------------------------------------------------------|--------|
| <b>SFC</b> method score:                                 | 50.29  |
| Optimised small scale <b>HPLC</b> method (not performed) | 202.10 |

*Assuming the same number of injections on both techniques*



# Conclusions

- Not surprisingly, SFC is typically « greener » for prep
- If you can push to go without co-solvent, even better (100% CO<sub>2</sub>)
- AMGS calculator is a useful tool to quickly compare different approaches in terms of sustainability





Thanks to the team!

Virginie Gonnord  
Aurélie Bich  
Noémie Viller

Thanks for your attention!

