

P/N 920-00138-03 Article Structure Description	
Syringe	SYRC HS5.0-23-5
Sample	Gas tight syringe for headspace injection
Maximum operation Temperature	110°C
Volume	5.0ml
Needle	Fixed needle
Needle length	51mm
Needle outer diameter	23 gauge (0,64mm)
Point style needle tip	Pointstyle 5 (conical needle tip and side hole)
Specials	Hole on top of syringe glass body for automatic gas flushing
Additional remarks	Temperature restriction due to glue used to fix the needle in the glass body. For high temperature applications same syringe with high temperature glue is available.
Packing	1 pcs