

P/N 980-04382 Article Structure Description	
Vial type,	1.5ml Screw Neck Vial ND8
Size	32 x 11.6mm
Glass material	Clear glass
Glass quality	1st hydrol. class
Volume	1.5 ml
Additional remarks	Vial with label and filling lines
Packing	1000 pcs



P/N 980-04383 Article Structure Description	
Cap	Black PP Thread Cap
Cap center hole size	5.5 mm
Septum type	Silicone white/PTFE red ultra clean
Septum thickness	1.3 mm
Septum hardness	45° shore A
Temperature resistant	-60 to 200°C
Additional remarks	
Packing	1000 pcs

